

Product datasheet for RC229607L3V

CXorf40A (EOLA1) (NM_001171907) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	EOLA1
Synonyms:	CXorf40; CXorf40A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001171907
ORF Size:	474 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC229607).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001171907.1
RefSeq Size:	2455 bp
RefSeq ORF:	477 bp
Locus ID:	91966
UniProt ID:	Q8TE69
Cytogenetics:	Xq28



MW: 17.9 kDa

Gene Summary: May have an important role of cell protection in inflammation reaction.[UniProtKB/Swiss-Prot Function]