

Product datasheet for **RC229613**

BEX2 (NM_001168399) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: BEX2 (NM_001168399) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: BEX2
Synonyms: BEX1; DJ79P11.1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC229613 representing NM_001168399
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCAGAAAATGGTGGTTTGCAGGGCCAAAGTGTGCGGCGACGCACCTCACGTCGAGAATCGGGAGGAGG
AGACTGCAAGGATAGGCCAGGAGTAATGGAGTCAAAGAGGAACGAGCGTTAAACAATCTCATCGTGGA
AAATGTCAACCAGGAAAATGATGAAAAGATGAAAAGGAGCAAGTTGCTAATAAAGGGGAGCCCTTGCC
CTACCTTTGAATGTTAGTGAATACTGTGTGCTAGAGGAAACCGTAGGCGGTTCCGCGTTAGGCAGCCCA
TCTGCAGTATAGATGGGACATAATGCATAGGCTTGGAGAGCCACAGGCAAGGATGAGAGAGGAGAATAT
GGAAAGGATTGGGGAGGAGGTGAGACAGCTGATGGAAAAGCTGAGGGAAAAGCAGTTGAGTCATAGTTTG
CGGGCAGTCAGCACTGATCCCCCTACCATGACCATCACGATGAGTTTTGCCTTATGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC229613 representing NM_001168399
Red=Cloning site Green=Tags(s)

MQKMVVC GAKCCGDAPHVENREEETARIGPVMESKEERALNNLIVENVNQENDEKDEKEQVANKGEPLA
LPLNVSEYCVPRGNRRFRVRQPILQYRWDIMHRLGEPQARMREENMERIGEEVRLMEKLRKQLSHSL
RAVSTDPHHDHDEFCLMP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites: SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_001168399

ORF Size: 480 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001168399.1, NP_001161871.1</u>
RefSeq ORF:	483 bp
Locus ID:	84707
UniProt ID:	<u>Q9BXY8</u>
Cytogenetics:	Xq22.2
MW:	19.2 kDa
Gene Summary:	This gene belongs to the brain expressed X-linked gene family. The encoded protein interacts with the transcription factor LIM domain only 2 in a DNA-binding complex that recognizes the E-box element and promotes transcription. This gene has been found to be a tumor suppressor that is silenced in human glioma. In breast cancer cells, this gene product modulates apoptosis in response to estrogen and tamoxifen, and enhances the anti-proliferative effect of tamoxifen. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009]