

Product datasheet for **RC229627L3V**

RRAS2 (NM_001177314) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RRAS2 (NM_001177314) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RRAS2
Synonyms:	NS12; TC21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001177314
ORF Size:	507 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229627).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001177314.1 , NP_001170785.1
RefSeq ORF:	510 bp
Locus ID:	22800
UniProt ID:	P62070
Cytogenetics:	11p15.2
Protein Families:	Druggable Genome
Protein Pathways:	MAPK signaling pathway, Regulation of actin cytoskeleton, Tight junction
MW:	20.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the R-Ras subfamily of Ras-like small GTPases. The encoded protein associates with the plasma membrane and may function as a signal transducer. This protein may play an important role in activating signal transduction pathways that control cell proliferation. Mutations in this gene are associated with the growth of certain tumors. Pseudogenes of this gene are found on chromosomes 1 and 2. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Apr 2010]