

## Product datasheet for **RC229635L4V**

### **C9orf25 (FAM219A) (NM\_001184942) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	C9orf25 (FAM219A) (NM_001184942) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C9orf25
Synonyms:	C9orf25
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001184942
ORF Size:	519 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229635).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001184942.1</a>
RefSeq ORF:	522 bp
Locus ID:	203259
UniProt ID:	<a href="#">Q8IW50</a>
Cytogenetics:	9p13.3
MW:	19.5 kDa
Gene Summary:	The protein encoded by this gene has homologs that have been identified in mouse, macaque, etc organisms. Multiple alternatively spliced transcript variants that encode different protein isoforms have been described for this gene. [provided by RefSeq, Dec 2010]


[View online »](#)