

## Product datasheet for **RC229645L3V**

### HDHD1A (PUDP) (NM\_001178136) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HDHD1A (PUDP) (NM_001178136) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HDHD1A
Synonyms:	DXF68S1E; FAM16AX; GS1; HDHD1; HDHD1A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001178136
ORF Size:	555 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229645).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001178136.1</a>
RefSeq ORF:	558 bp
Locus ID:	8226
UniProt ID:	<a href="#">Q08623</a>
Cytogenetics:	Xp22.31
MW:	21 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the haloacid dehalogenase-like (HAD) hydrolase superfamily. The encoded protein has no known biological function. This gene has a pseudogene on chromosome 1. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2010]