

Aquaporin 1 (AQP1) (NM_001185060) Human Tagged Lenti ORF Clone


Product Type:	Expression Plasmids
Product Name:	Aquaporin 1 (AQP1) (NM_001185060) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Aquaporin 1
Synonyms:	AQP-CHIP; CHIP28; CO
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229650).
Restriction Sites:	Sgfl-Mlul
Cloning Scheme:	

SgfI ORF *MluI*
 ... GCG ATC GC C ATG --- // --- NNN⁺ ACG CGT ...

[illegible]

* The last codon before the Stop codon of the ORE.

[View online »](#)

Plasmid Map:


ACCN: NM_001185060

ORF Size: 558 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001185060.1](#), [NP_001171989.1](#)

RefSeq ORF: 561 bp

Locus ID: 358

Cytogenetics:	7p14.3
Protein Families:	Druggable Genome, Ion Channels: Other, Transmembrane
MW:	20.4 kDa
Gene Summary:	This gene encodes a small integral membrane protein with six bilayer spanning domains that functions as a water channel protein. This protein permits passive transport of water along an osmotic gradient. This gene is a possible candidate for disorders involving imbalance in ocular fluid movement. [provided by RefSeq, Aug 2016]