

## Product datasheet for **RC229656L4V**

### CD99L2 (NM\_001184808) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CD99L2 (NM_001184808) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CD99L2
Synonyms:	CD99B; MIC2L1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001184808
ORF Size:	567 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229656).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001184808.1</a> , <a href="#">NP_001171737.1</a>
RefSeq ORF:	570 bp
Locus ID:	83692
UniProt ID:	<a href="#">Q8TCZ2</a>
Cytogenetics:	Xq28
Protein Families:	Transmembrane
MW:	20.5 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a cell-surface protein that is similar to CD99. A similar protein in mouse functions as an adhesion molecule during leukocyte extravasation. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]