

## Product datasheet for **RC229733**

### CLM 9 (CD300LG) (NM\_001168324) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CLM 9 (CD300LG) (NM_001168324) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CLM 9
Synonyms:	CLM-9; CLM9; NEPMUCIN; TREM-4; TREM4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC229733 representing NM_001168324 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**

ATGCGGCTTCTGGTCTGCTATGGGGTTCCTGCTGCTCCAGGTTATGAAGCCCTGGAGGGCCAGAGGAAATCAGCGGGTTCGAAGGGGACACTGTGTCCCTGCAGTGCACCTACAGGGAAGAGCTGAGGGACCACCGAAAGTACTGGTGCAGGAAGGGTGGGATCCTTCTCTCGCTGCTCTGGCACCATCTATGCAGAAGAAGAAAGCCAGGAGACAATGAAGGGCAGGGTGTCCATCCGTGACAGCCGCCAGGAGCTCTCGCTCATTGTGACCCGTGTGAACTCACCTGCAAGACGCTGGGGAGTACTGGTGTGGGGTCGAAAACGGGGCCCGATGAGCTTTACTGATCTCTGTTTCAGGGAGCTCCCGCCCCCATGCAGCTGGACTCCACCTCAGCAGGACACCAGTCCAGCTCTCAGCAGTGGCAGCTCTAAGCCCAGGGTGTCCATCCCGATGGTCCGCATATGGCCCCAGTCCCTGGTGTGCTGAGCCTTCTGTGAGCCGAGGCCTGATCGCCTTCTGCAGCCACCTGCTCCTGTGGAGAAAGGAAGCTCAACAGGCCACGGAGACACAGAGGAACGAGAAGTTCTGCCTCTCACGCTTGAAATCACAGAACATCAGAGCTGGATCCAGATGAAGCCTTTGAGATCTTAGCCTGACTTTTCTCTCAGTT

**ACGCGT**ACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC229733 representing NM\_001168324  
Red=Cloning site Green=Tags(s)

MRLLVLLWGCLLLPGYEALEGPPEISGFEGDVTSLQCTYREELRDHRKYWCRKGGILFSRCSGTIYAE  
 EQETMKGRVSIIRDSRQELSLIVTLWNLTLQDAGEYWCVEKRPDESLLISLFFVPGSSRPPMQLDSTSA  
 EDTSPALSSGSSKPRVSIIPMVRILAPVLLSLLSAAGLIAFCSHLLLWRKEAQQATETQRNEKFLSRL  
 NHRTSELDPDEAFEISSLTFLSV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_001168324

**ORF Size:** 699 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001168324.2](#)

**RefSeq ORF:** 702 bp

**Locus ID:** 146894

**UniProt ID:** [Q6UXG3](#)


**Cytogenetics:** 17q21.31

**Protein Families:** Druggable Genome, Transmembrane

**MW:** 26.5 kDa

**Gene Summary:** Members of the CD300 (see MIM 606786)-like (CD300L) family, such as CD300LG, are widely expressed on hematopoietic cells. All CD300L proteins are type I cell surface glycoproteins that contain a single immunoglobulin (Ig) V-like domain (Takatsu et al., 2006 [PubMed 16876123]).[supplied by OMIM, Mar 2008]

**Product images:**



Circular map for RC229733