

## Product datasheet for **RC229757L4V**

### PEX11B (NM\_001184795) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PEX11B (NM_001184795) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PEX11B
Synonyms:	PEX11-BETA; PEX14B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001184795
ORF Size:	735 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229757).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001184795.1</a> , <a href="#">NP_001171724.1</a>
RefSeq ORF:	738 bp
Locus ID:	8799
UniProt ID:	<a href="#">O96011</a>
Cytogenetics:	1q21.1
Protein Families:	Transmembrane
MW:	27.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene facilitates peroxisomal proliferation and interacts with PEX19. The encoded protein is found in the peroxisomal membrane. Several transcript variants, some protein-coding and some not protein-coding, have been found for this gene. [provided by RefSeq, Dec 2012]