

Product datasheet for **RC229768**

OCIAD1 (NM_001168254) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OCIAD1 (NM_001168254) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OCIAD1
Synonyms:	ASRIJ; OCIA; TPA018
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC229768 representing NM_001168254 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTGC GGAAAGATGAATGGGAGGGCTGATTTTCGAGAGCCGAATGCAGAGGTTCCAAGACCAATTC
CCCACATAGGGCCTGATTACATTC AACAGAGGAAGAAAGGAGAGTCTTCGAGAATGCAATGATGAAAG
CTTCTGGTTCAGATCTGTGCCTTTGGCTGCAACAAGTATGTTGATTACTCAAGGATTAATTAGTAAAGGA
ATACTTTCAAGTCATCCCAAATATGGTTCCATCCCTAACTTATACTTGCTTGATCATGGGATACTTTG
CTGAAAACTTTCTTATGTGAAAACTTGCCAAGAGAAATTC AAGAACTTGAAAAATCCCCCCTTGAGAG
AGCTTTACGATCAGGACAAGCAGCAGATCTTACCACCTGGGCACTATTATCAAAAGTCAAATATGAC
TCAAGTGTGAGTGGTCAATCATCTTTTGTGACATCCCCAGCAGCAGACAACATAGAAATGCTTCCTCATT
ATGAGCCAATTCATTCACTTCTTATGAATGAATCTGCTCCCACTGGTATTACTGATCATATTGTCCA
AGGACCTGATCCCAACCTTGAAGAAAGTCTAAAAGAAAAAATATTACATATGAGGAATTAAGGAATAAG
AACAGAGATCATATGAAGTATCTTTAACACAAAAGACTGACCCCTCAGTCAGGCCTATGCATGAAAGAG
TGCCAAAAAAGAAGTCAAAGTAAACAAGTATGGAGATACTTGGGATGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC229768 representing NM_001168254
 Red=Cloning site Green=Tags(s)

MECGKMNGRADFREPNAEVPRPIPHIGPDYIPTEEERRVFAECNDESFWFRSVPLAATSMITQGLISKG
 ILSSHPKYGSIPKLILACIMGYFAGKLSYVKTCQEKFKKLENSPLGEALRSGQARRSSPPGHYYQSKYD
 SSVSGQSSSVTSPAADNIEMLPHYEPIPFSSSMNESAPTGITDHIIVQGPDPNLEESPKRKNITYEELRNK
 NRESYEVSLTQKTDPVSRPMPHERVPKKEVKVNYGDTWDE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001168254

ORF Size: 750 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001168254.1](#), [NP_001161726.1](#)

RefSeq ORF: 753 bp


Locus ID: 54940

Cytogenetics: 4p11

MW: 28.6 kDa

Gene Summary: Maintains stem cell potency (By similarity). Increases STAT3 phosphorylation and controls ERK phosphorylation (By similarity). May act as a scaffold, increasing STAT3 recruitment onto endosomes (By similarity). Involved in integrin-mediated cancer cell adhesion and colony formation in ovarian cancer (PubMed:20515946).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC229768