

## Product datasheet for **RC229856L4V**

### PEX2 (NM\_001172086) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PEX2 (NM_001172086) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PEX2
Synonyms:	PAF1; PBD5A; PBD5B; PMP3; PMP35; PXMP3; RNF72; ZWS3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001172086
ORF Size:	915 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229856).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001172086.1</a> , <a href="#">NP_001165557.1</a>
RefSeq ORF:	918 bp
Locus ID:	5828
UniProt ID:	<a href="#">P28328</a>
Cytogenetics:	8q21.13
Protein Families:	Druggable Genome, Transmembrane
MW:	35.3 kDa



[View online »](#)

**Gene Summary:**

This gene encodes an integral peroxisomal membrane protein required for peroxisome biogenesis. The protein is thought to be involved in peroxisomal matrix protein import. Mutations in this gene result in one form of Zellweger syndrome and infantile Refsum disease. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Jul 2008]