

## Product datasheet for **RC229876L4V**

### G protein alpha Inhibitor 2 (GNAI2) (NM\_001166425) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	G protein alpha Inhibitor 2 (GNAI2) (NM_001166425) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GNAI2
Synonyms:	GIP; GNAI2B; H_LUCA15.1; H_LUCA16.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001166425
ORF Size:	954 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229876).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001166425.1</a>
RefSeq ORF:	957 bp
Locus ID:	2771
UniProt ID:	<a href="#">P04899</a>
Cytogenetics:	3p21.31
Protein Families:	Druggable Genome



[View online »](#)

<b>Protein Pathways:</b>	Axon guidance, Chemokine signaling pathway, Gap junction, Leukocyte transendothelial migration, Long-term depression, Melanogenesis, Progesterone-mediated oocyte maturation, Tight junction
<b>MW:</b>	36.9 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is an alpha subunit of guanine nucleotide binding proteins (G proteins). The encoded protein contains the guanine nucleotide binding site and is involved in the hormonal regulation of adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2013]