

Product datasheet for RC229912

PPAR delta (PPARD) (NM_001171820) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPAR delta (PPARD) (NM_001171820) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PPAR delta
Synonyms:	FAAR; NR1C2; NUC1; NUCI; NUCII; PPARB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC229912 representing NM_001171820 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGCAGCCACAGGAGGAAGCCCTGAGGTCCGGGAAGAGGAGGAGAAAGAGGAAGTGGCAGAGGCAG
AAGGAGCCCCAGAGCTCAATGGGGGACCACAGCATGCACCTTCTCCAGCAGCTACACAGCTATCCGTTT
TGGTCCGATGCCGGAGGCTGAGAAGAGGAAGCTGGTGGCAGGGCTGACTGCAAACGAGGGGAGCCAGTAC
AACCCACAGGTGGCCGACCTGAAGGCCTTCTCCAAGCACATCTACAATGCCTACCTGAAAACTTCAACA
TGACAAAAAGAAGGCCCGCAGCATCTCACCGCAAAGCCAGCCACACGCGCCCTTTGTGATCCACGA
CATCGAGACATTGTGGCAGGCAGAGAAGGGGCTGGTGTGGAAGCAGTTGGTGAATGGCCTGCCTCCCTAC
AAGGAGATCAGCGTGCACGTCTTCTACCGCTGCCAGTGCACCACAGTGGAGACCGTGCGGGAGCTCACTG
AGTTCCGCAAGAGCATCCCCAGCTTCCAGCAGCCTTCTCTCAACGACCAGGTTACCTTCTCAAGTATGG
CGTGCACGAGGCCATCTTCGCCATGCTGGCCTCTATCGTCAACAAGGACGGGCTGCTGGTAGCCAACGGC
AGTGGCTTTGTCAACCGTGAGTTCCTGCGCAGCCTCCGCAAACTTTCAGTGATATCATTGAGCCTAAGT
TTGAATTTGCTGTCAAGTTCAACGCCCTGGAACCTGATGACAGTGACCTGGCCCTATTATTGCGGCCAT
CATTCTGTGTGGAGACCGCCAGGCCTCATGAACGTTCCACGGGTGGAGGCTATCCAGGACACCATCCTG
CGTGCCCTCGAATTCACCTGCAGGCCAACCCCTGATGCCAGTACCTTTCCCAAGCTGCTGCAGA
AGATGGCTGACCTGCGGCAACTGGTCACCGAGCACGCCAGATGATGCAGCGGATCAAGAAGACCCGAAAC
CGAGACCTCGCTGCACCCTCTGCTCCAGGAGATCTACAAGGACATGTAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC229912 representing NM_001171820
Red=Cloning site Green=Tags(s)

MEQPQEEAPEVREEEKEEVAEAEAGAPELNGGPQHLPSSSYTAIRFGRMPEAEKRKLVAGLTANEGSQY
 NPQVADLKAFSKHIYNAYLKNFNMTKKKARSILTGKASHTAPFVIHDIETLWQAEKGLVWKLVNGLPPY
 KEISVHVFYRCQCTTVETVRELTEFAKSIPSFSSLFLNDQVTLKYGVEAIFAMLASIVNKDGLLVANG
 SGFVTREFLRSLRKPFSDIIEPKFEFAVKFNALELDDSDLALFIAAIILCGDRPGLMNVPRVEAIQDTIL
 RALEFHLQANHPDAQYLFPKLLQKMADLRQLVTEHAQMMQRIKKTETETSLHPLLQEIKDMY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8053_e10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001171820

ORF Size: 1029 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001171820.1](#), [NP_001165291.1](#)

RefSeq ORF: 1032 bp

Locus ID: 5467

UniProt ID: [Q03181](#)

Cytogenetics: 6p21.31


Protein Families: Druggable Genome, Nuclear Hormone Receptor, Transcription Factors

Protein Pathways: Acute myeloid leukemia, Pathways in cancer, PPAR signaling pathway, Wnt signaling pathway


MW: 39.3 kDa

Gene Summary: This gene encodes a member of the peroxisome proliferator-activated receptor (PPAR) family. The encoded protein is thought to function as an integrator of transcriptional repression and nuclear receptor signaling. It may inhibit the ligand-induced transcriptional activity of peroxisome proliferator activated receptors alpha and gamma, though evidence for this effect is inconsistent. Expression of this gene in colorectal cancer cells may be variable but is typically relatively low. Knockout studies in mice suggested a role for this protein in myelination of the corpus callosum, lipid metabolism, differentiation, and epidermal cell proliferation. Alternative splicing results in multiple transcript variants encoding distinct protein isoforms. [provided by RefSeq, Aug 2017]

Product images:



Circular map for RC229912



Western blot validation of overexpression lysate (Cat# [LY432912]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC229912 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).