

Product datasheet for **RC229923A1V**

Human PP11 (ENDOU) (NM_001172440) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PP11 (ENDOU) (NM_001172440) AAV Particle
Tag:	Myc-DDK
Symbol:	PP11
Synonyms:	P11; PP11; PRSS26
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC229923 representing NM_001172440 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**GCGATCGCC**

ATGAGGGCCTGCATCTCCCTGGTATTGGCCGTGCTGTGTGGCCTGGCCTGGGCTGACTTGTACTCGGCACCCACCTCTGCCAGGGCCGCTGCTACGAAGCCTTTGACAAGCACCACCAATGTCCTGCAATGCCCGCTGCCAAGAGTTTGGGAACGCTGCAAGGATTTTGGAGCCTGTGTAGTGACCACGAGGTCTCCACAGCAGTGATGCCATAACAAAAGAGGAGATTGAGCAGCATCTCTGAGAAGATCTACAGGGCAGACACCAACAAAGCCAGAAAGGAAGACATCGTTCTCAATAGCCAAAACGTCATCTCCCCGTCAGAGACCAGAAACCAAGTGGATCGCTGCCAAAGCCACTTTCACCTTATGTCAATGAGAAGCTGTTCTCCAAGCCACCTATGCAGCCTTCATCAACCTCTCAACAACCTACCAGCGGGCAACAGGCCATGGGGAGCACTTCAGTGCCAGGAGCTGGCCGAGCAGGACGCTTCCCTCAGAGAGATCATGAAGACAGCAGTCATGAAGGAGCTTACAGCTTCCATCACCAGAATCGCTATGGCTCAGAGCAAGAGTTTGTGCGATGACTTGAAGAACATGTGGTTTGGGCTCTATTCAAGAGGCAATGAAGAGGGGGACTCGAGTGGCTTTGAACATGTCTTCAGGTGAGGTAAAAAAGGCAAGGTTACTGGCTCCATAACTGGATCCGCTTCTACCTGGAGGAGAAGGAGGTCTGGTTGACTATTACAGTCACATCTACGATGGGCCTTGGGATTCTACCCCGATGTGCTGGCAATGCAGTTCAACTGGGACGGCTACTATAAGGAAGTGGGCTCTGCTTTCATCGGCAGCAGCCCTGAGTTTGTAGTTTGCCTTACTCCCTGTGCTTCATCGCCAGGCCAGGCAAGTGTGCCAGTTAAGCCTGGGAGGATATCCCTTAGCTGTCCGGACATATACCTGGGA CAAGTCCACCTATGGGAATGGCAAGAAGTACATCGCCACAGCCTACATAGTGTCTTCCACC

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC229923 representing NM_001172440 Red=Cloning site Green=Tags(s) MRACISLVLAVLCGLAWADLYSAPTSCQGRCYEAFDKHHQCHCNARCQEFGNCKDFESLCSDFEVSHTSS DAITKEEIQSISEKIYRADTNKAQKEDIVLNSQNCISPSETRNQVDRCPKPLFTYVNEKLF SKPTYAAFI NLLNNYQRATGHGEHFS AQELAEQDAFLREIMKTAVMKELYSFLHHQNRYGSEQEFVDDLKMMWFLYSR GNEEGDSSGFEHVFSGEVKKGKVTGFHNWIRFYLEEKEGLVDYYSHIYDGPWDSYPDVLAMQFNWDGYK EVGSAFIGSSPEFEFALYSLCFIARPGKVCQLSLGGYPLAVRITYWTKSTYNGKKYIATAYIVSST TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001172440
ORF Size:	1041 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001172440.1
RefSeq ORF:	1044 bp
Locus ID:	8909
UniProt ID:	P21128
Cytogenetics:	12q13.11
MW:	40.1 kDa