

## Product datasheet for **RC229934A1V**

### Human PDHA1 (NM\_001173456) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human PDHA1 (NM\_001173456) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** PDHA1  
**Synonyms:** PDHA; PDHAD; PDHCE1A; PHE1A  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC229934 representing NM\_001173456  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGGAAGATGCTCGCCGCGTCTCCCGCGTGTCTGGCGTTCTCAGAAGCCGGCAAGCAGAGTGC  
TGGTAGCATCCCGTAATTTGCAAAATGATGCTACATTTGAAATTAAGAAATGTGACCTTCACCGGCTGGA  
AGAAGGCCCTCCTGTCAACAAGTGTCAACAGGAGGATGGGCTCAAATACTACAGGATGATGCAGAGT  
GTACGCCGAATGGAGTTGAAAGCAGATCAGCTGTATAAACAGAAAATTATTCGTGGTTTTCTGCACTTGT  
GTGATGGTCAGGAAGCTTGTGTGTGGCCTGGAGGCCGCATCAACCCACAGACCATCTCACAGC  
CTACCGGGCTCACGGCTTTACTTTACCCGGGGCCTTTCCGTCGAGAAAATTCTCGCAGAGCTTACAGGA  
CGAAAAGGAGGTTGTGCTAAAGGAAAAGGAGGATCGATGCACATGTATGCCAAGAATTCTACGGGGCA  
ATGGCATCGTGGGAGCGCAGGGCCAGATATTCGAAGCTTACAACATGGCAGCTTTGTGGAATTACCTTG  
TATTTTCATCTGTGAGAATAATCGCTATGGAATGGGAACGTCTGTTGAGAGAGCGGCAGCCAGCACTGAT  
TACTACAAGAGAGGCGATTTTCTTCTGGGCTGAGAGTGGATGGAATGGATATCCTGTGCGTCCGAGAGG  
CAACAAGGTTTGTGCTGCCTATTGTAGATCTGGGAAGGGCCCATCCTGATGGAGCTGCAGACTTACCG  
TTACCACGGACACAGTATGAGTGACCCTGGAGTCAGTTACCGTACACGAGAAGAAATTCAGGAAGTAAGA  
AGTAAGAGTGACCCTATTATGCTTCTCAAGGACAGGATGGTGAACAGCAATCTTGCCAGTGTGGAAGAAC  
TAAAGGAAATGATGTGGAAGTGAAGGAGATTGAGGATGCTGCCAGTTTGGCACGGCCGATCCTGA  
GCCACCTTTGGAAGAGCTGGGCTACCACATCTACTCCAGCGACCCACCTTTTGAAGTTCTGTGGTCCAAT  
CAGTGGATCAAGTTAAGTCAGTCAGT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Protein Sequence:</b>	>RC229934 representing NM_001173456 Red=Cloning site Green=Tags(s)  MRKMLAAVSRVLSGASQKPAASRVLVASRNFANDATFEIKKCDLHRLLEEGPPVTTVL TREDGLKYYRMMQT VRRMELKADQLYKQKIIRGFCHLCDGQEACCVGLEAGINPTDHLITAYRAHGFTFTRGLSVREILAEELTG RKGGCAKGGGSMHMYAKNFYGGNGIVGAQGQIFEAYNMAALWKLPCIF ICENNRVYGMGT SVERAAASTD YYKRGDFIPGLRVDGMDILCVREATRFAAAYCRSGKGPILMELQTYRYHGHSMSPGVS YRTREEIQEVR SKSDPIMLLKDRMVNSNLA SVEELKEIDVEVRKEIEDAAQFATADPEPPLEELGYHIYSSDPPFEVRGAN QWIKFKSVS  TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
<b>Species:</b>	Human
<b>Serotype:</b>	AAV-2
<b>ACCN:</b>	NM_001173456
<b>ORF Size:</b>	1077 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<a href="#">NM_001173456.1</a>
<b>RefSeq ORF:</b>	1080 bp
<b>Locus ID:</b>	5160
<b>UniProt ID:</b>	<a href="#">P08559</a>
<b>Cytogenetics:</b>	Xp22.12
<b>MW:</b>	40.6 kDa