

Product datasheet for **RC229950L4V**

OGR1 (GPR68) (NM_001177676) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	OGR1 (GPR68) (NM_001177676) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OGR1
Synonyms:	A12A6; GPR12A; OGR1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001177676
ORF Size:	1095 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229950).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001177676.1
RefSeq ORF:	1098 bp
Locus ID:	8111
UniProt ID:	Q15743
Cytogenetics:	14q32.11
Protein Families:	Druggable Genome, GPCR, Transmembrane
MW:	41.5 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a G protein-coupled receptor for sphingosylphosphorylcholine. The encoded protein is a proton-sensing receptor, inactive at pH 7.8 but active at pH 6.8. Mutations in this gene are a cause of amelogenesis imperfecta. [provided by RefSeq, Feb 2017]