

Product datasheet for **RC229966L4V**

SRPX (NM_001170752) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SRPX (NM_001170752) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SRPX
Synonyms:	DRS; ETX1; HEL-S-83p; SRPX1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001170752
ORF Size:	1137 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC229966).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001170752.1 , NP_001164223.1
RefSeq ORF:	1140 bp
Locus ID:	8406
UniProt ID:	P78539
Cytogenetics:	Xp11.4
MW:	42.2 kDa
Gene Summary:	May be involved in phagocytosis during disk shedding, cell adhesion to cells other than the pigment epithelium or signal transduction.[UniProtKB/Swiss-Prot Function]



[View online »](#)