

Product datasheet for **RC229969A1V**

Human FCRLM1 (FCRLA) (NM_001184866) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FCRLM1 (FCRLA) (NM_001184866) AAV Particle
Tag:	Myc-DDK
Symbol:	FCRLM1
Synonyms:	FCRL; FCRL1; FCRLb; FCRLc1; FCRLc2; FCRLd; FCRLe; FCRLM1; FCRLX; FCRX; FRED
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC229969 representing NM_001184866 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTTGAAGAAATCAGTGTGGGGTTGCAGGAGACCTAAACACAGTCACCATGAAGCTGGGCTGTGTCC
 TCATGGCCTGGGCCCTACCTTTCCCTTGGTGTGCTCTGGGTGGCCAGATGCTACTGGCAGCTGGATG
 TCATGCCGCTGCCAGTTTGGAGACGCTGCAGTGTGAGGGACCTGTCTGCACTGAGGAGAGCAGCTGCCAC
 ACGGAGGATGACTTGACTGATGCAAGGGAAGCTGGCTTCCAGGTCAAGGCCTACACTTTCAGTGAACCT
 TCCACCTGATTGTGCTATGACTGGCTGATCCTCAAGGTCCAGCCAAGCCAGTTTTTGAAGGGGACCT
 GCTGGTTCTGCGCTGCCAGGCCTGGCAAGACTGGCCACTGACTCAGGTGACCTTCTACCGAGATGGCTCA
 GCTCTGGGTCCCCCGGGCCTAACAGGGAATTCTCCATCACCGTGGTACAAAAGGCAGACAGCGGGCACT
 ACCACTGCAGTGGCATCTTCCAGAGCCCTGGTCTGGGATCCCAGAAACAGCATCTGTTGTGGCTATCAC
 AGTCCAAGAACTGTTTCCAGCGCAATTCTCAGAGCTGTACCCTCAGCTGAACCCCAAGCAGGAAGCCCC
 ATGACCCTGAGTTGTGACACAAAGTTGCCCTGCAGAGGTGAGTGGCCGCTCCTCTTCTCTTCTACA
 AGGATGGAAGGATAGTGCAAAGCAGGGGGCTCTCCTCAGAATTCCAGATCCCCACAGCTTCAGAAGATCA
 CTCCGGGTCACTGTTGTGAGGCAGCCACTGAGGACAACCAAGTTTGGAAACAGAGCCCCCAGCTAGAG
 ATCAGAGTGCAGGTGCTTCCAGCTCTGCTGCACCTCCACATTGAATCCAGCTCCTCAGAAATCAGCTG
 CTCCAGGAAGTCTCCTGAGGAGGCCCTGGGCTCTGCTCCGCCGCCAACCCATCTTCTGAGGATCC
 AGGCTTTTCTTCTCTGTTGATGCCAGATCCTCATCTGTATACCAGATGGGCTTCTTCTCAAACAC
 ATGCAGGATGTGAGAGTCTCCTCGGTACCTGCTCATGGAGTTGAGGGAATTATCTGGCCACCGGAAGC
 CTGGGACCACAAAGGCTACTGCTGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence:	>RC229969 representing NM_001184866 Red=Cloning site Green=Tags(s) MLKKISVGVDLNTVTMKLGCVLMAWALYLSLGVLWVAQMLLAAGCHAAASFETLQCEGPVCTEESCH TEDDLTDAREAGFQVKAYTFSEPFHLIVSYDWLILQGPAKPVFEGDLLVLRQAWQDWPLTQVTFYRDGS ALGPPGPNREFSITVVQKADSGHYHCSGIFQSPGPGIPETASVVAITVQELFPAPILRAVPSAEPQAGSP MTLSCQTKLPLQRSAARLLFSFYKDGRIVQSRGLSSEFQIPTASEDHSGSYWCEAATEDNQVWKQSPQLE IRVQGASSAAPPPTLNPAPQKSAAPGTAPEEAPGPLPPPTPSSDPGFSSPLGMPDPLHYHQMGLLLKH MQDVRVLLGHLLMELRELSGHRKPGTTKATAE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001184866
ORF Size:	1146 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001184866.1</u> , <u>NP_001171795.1</u>
RefSeq ORF:	1098 bp
Locus ID:	84824
UniProt ID:	<u>Q7L513</u>
Cytogenetics:	1q23.3
MW:	41.6 kDa