

## Product datasheet for **RC229980L3V**

### **FAM46D (TENT5D) (NM\_001170574) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

|                                  |                                      |
|----------------------------------|--------------------------------------|
| <b>Product Type:</b>             | Lentiviral Particles                 |
| <b>Symbol:</b>                   | TENT5D                               |
| <b>Synonyms:</b>                 | CT1.26; CT112; FAM46D                |
| <b>Mammalian Cell Selection:</b> | Puromycin                            |
| <b>Vector:</b>                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| <b>Tag:</b>                      | Myc-DDK                              |
| <b>ACCN:</b>                     | NM_001170574                         |
| <b>ORF Size:</b>                 | 1167 bp                              |

**ORF Nucleotide Sequence:** The ORF insert of this clone is exactly the same as(RC229980).

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**RefSeq:** [NM\\_001170574.1](#), [NP\\_001164045.1](#)

**RefSeq ORF:** 1170 bp

**Locus ID:** 169966

**UniProt ID:** [Q8NEK8](#)

**Cytogenetics:** Xq21.1

**MW:** 44.9 kDa



**Gene Summary:**

Antibodies against the protein encoded by this gene were found only in plasma from cancer patients. While it may be a target for immunotherapy, the function of this gene is unknown.  
[provided by RefSeq, Dec 2009]