

Product datasheet for **RC230031L4V**

MEPE (NM_001184695) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | MEPE (NM_001184695) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | MEPE |
| Synonyms: | OF45 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001184695 |
| ORF Size: | 1236 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC230031). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001184695.1 |
| RefSeq Size: | 2113 bp |
| RefSeq ORF: | 1239 bp |
| Locus ID: | 56955 |
| UniProt ID: | Q9NQ76 |
| Cytogenetics: | 4q22.1 |
| Protein Families: | Secreted Protein |
| MW: | 45.3 kDa |



[View online »](#)

Gene Summary:

This gene encodes a secreted calcium-binding phosphoprotein that belongs to the small integrin-binding ligand, N-linked glycoprotein (SIBLING) family of proteins. Members of this family are components of the extracellular matrix of bone and dentin and regulate bone mineralization. Deficiency of a similar protein in mouse results in increased bone mass. Mice lacking this gene are resistant to aging-related trabecular bone loss. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2014]