

Product datasheet for **RC230036L3V**

ZNF641 (NM_001172682) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ZNF641
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001172682
ORF Size:	1245 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230036).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001172682.1
RefSeq ORF:	1248 bp
Locus ID:	121274
UniProt ID:	Q96N77
Cytogenetics:	12q13.11
Protein Families:	Transcription Factors
MW:	47 kDa



Gene Summary:

Transcriptional activator. Activates transcriptional activities of SRE and AP-1.[UniProtKB/Swiss-Prot Function]