

Product datasheet for **RC230056L3V**

VSIG1 (NM_001170553) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | VSIG1 (NM_001170553) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | VSIG1 |
| Synonyms: | 1700062D20Rik; dj889N15.1; GPA34 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001170553 |
| ORF Size: | 1269 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC230056). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001170553.1 , NP_001164024.1 |
| RefSeq ORF: | 1272 bp |
| Locus ID: | 340547 |
| UniProt ID: | Q86XK7 |
| Cytogenetics: | Xq22.3 |
| Protein Families: | Transmembrane |
| MW: | 46.3 kDa |



[View online »](#)

Gene Summary:

This gene encodes a member of the junctional adhesion molecule (JAM) family. The encoded protein contains multiple glycosylation sites at the N-terminal region, and multiple phosphorylation sites and glutamic acid/proline (EP) repeats at the C-terminal region. The gene is expressed in normal stomach and testis, as well as in gastric, esophageal and ovarian cancers. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009]