

Product datasheet for **RC230057**

ZNF641 (NM_001172681) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF641 (NM_001172681) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF641
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC230057 representing NM_001172681 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTTTTCAAGACAGACAGCAGCGCTGGGGACAGGATGGGAATCAATGAATGTACAGCTGGATGGAGCAG
 AGCCCCAGGTGGAAGGGGAAGCCAGGAAGAGCGGCCATGGAGAACAGTACCAGGACCTCTGGAGCACCT
 GTGCTGTGACCTTGAAGAGGAGCCACAGTCCCTTCAGGAGAAGGCTCAGTCTGCTCCCTGGGTTCTTGCA
 ATCCCCAGGAGGGGAACACTGGAGACTGGGAGATGGCAGCTGCATTCTTGGCGCTGGATCACAGGGCC
 TGGTAACCATCAAGGATGTGCTACTGTGCTTCTCTCAGGAGGAGTGGCGGAGCCTGGACCCCTCTCAGAC
 AGACTTTTATGGAGAATATGTCATGCAGGAAAACCTGTGGGATAGTAGTCTCTCTGAGATTTCCAATTCCC
 AAAGTGGACATGCTTTCTCAACTAGAAGGAGGAGAAGAACAATGGGTCCCTGACCCCCAGGACTTAGAGG
 AGAGGGACATTCTGAGGGTCACATATACAGGAGATGGAAGTGAACATGAGGGGGATACCCCTGAAGTAGA
 AGCAGAACCTCCAGAAATGTTATCCAGCGTGTCTGAAGATACTGTTCTCTGGAACCCGGAGCATGATGAG
 AGCTGGGATTCCATGCCAGCAGCTCCAGAGGAATGCTCCTGGGGCCCCCTTTTCTTCAGGAAGATAGTT
 TCTCAACCTGCTGTGTAGCACAGAGATGGATTCCCTGTTAAGACCCACACATGCCCCAGTGTGGGAA
 ACAGTTTGTATGGGGTTCCACCTTGCCAGGCATCAACAAACACACACTGGGAGAGACCCTACAGCTGC
 CTAAGTGTGAGAAGACCTTTGGGCGAAGACATCACCTCATCAGGCACCAGAAAACCCACCTACATGACA
 AGACCAGCAGGTGCTCTGAGTGTGGTAAGAATTTCCGATGCAACTCCCATCTGGCCAGCCACCAGAGAGT
 GCATGCAGAAGGCAATCCTGCAAAGGCCAAGAGGTTGGAGAGAGCCCTGGCACAAAGAAACGGCAGCGT
 GCCCCACCAAGTGCCAAAGTGTACGTGTGCACTGAATGTGGGAAGAGCTTTGGCCGAAGGCACCACCTTG
 TGAGACTGGCTGACCCACACTGGGGAGAAGCCCTTCCAGTGCCTCGCTGTGAGAAGAGCTTTGGCCG
 AAAACATCACCTGGACAGGCACCTGCTCACCCACCAGGACAAAGTCCCCGGAACAGCTGGGACAGAGGA
 ACATCTGTCTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC230057 representing NM_001172681
 Red=Cloning site Green=Tags(s)


MLSEQTAALGTGWESMNVLQDGAEPQVERGSQEERPWRTVPGPLEHLCCDLEEPQSLQEKAQSAPWVPA
 IPQEGNTGDWEMAAALLAAGSQGLVTIKDVSLCFSQEEWRSLDPSQTDIFYGEYVMQENCIVVSLRFP
 KLDMLSQLGEGEEQWVPDPQDLEERDILRVITYTGDGSEHEGDTPELEAEPRLMLSSVSEDTVLWNPEHDE
 SWDSMPSSSRGMLLGPPFLQEDSFNLLCSTEMDSLRLPHTCPQCGKQFVWGSHLARHQQTHTGERPYSC
 LKCEKTFGRRHHLIRHQKTHLHDKTSRCSECGKNFRCNHSLASHQVRVHAEGKSCKGQEVGESPGTRKRQR
 APPVPKCHVCTECGKSFGRHHLVRHWLTHTEKPFQCPRCEKSFGRKHHLDRHLLTHQGQSPRNSWDRG
 TSVF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001172681

ORF Size: 1272 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001172681.2

RefSeq Size: 4371 bp

RefSeq ORF: 1275 bp

Locus ID: 121274

UniProt ID: Q96N77


Cytogenetics: 12q13.11

Protein Families: Transcription Factors

MW: 48 kDa

Gene Summary: Transcriptional activator. Activates transcriptional activities of SRE and AP-1. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC230057