

Product datasheet for **RC230067**

ARL13B (NM_001174150) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARL13B (NM_001174150) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ARL13B
Synonyms:	ARL2L1; JBTS8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC230067 representing NM_001174150
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTTTCAGTCTGATGGCCAGTTGCTGCGGCTGGTTCAAGCGGTGGCGGGAGCCTGTCAGAAAGGTGACTC
 TTTTGATGGTGGGACTTGATAATGCTGGTAAAACCGCAACAGCAAGGGAATCCAAGGAGAATACCTGA
 AGATGTAGCTCCTACTGTTGGATTTTCAAAAATTAACCTTAGACAAGGAAAGTTGAAGTACCATCTTT
 GACTTGGGAGGTGGAATAAGAATTCGGGAATCTGGAAGAATTACTATGCTGAATCCTATGGGGTAATAT
 TTGTTGTGGATTCAGTGATGAAGAGAGAATGGAAGAGACAAAAGAGGCTATGTCAGAAATGCTAAGACA
 TCCTAGGATATCGGGAAGCCTATATTGGTGTGGCAAATAACAAGATAAAGAAGGAGCTTTAGGAGAA
 GCTGATGTCATTGAATGTCTATCTCTGGAAAAATTGGTCAATGAGCACAAGTGCCTGTGTCAGATAGAAC
 CATGTTACAGCAATCTCGGGTATGGAAAGAAAATTGACAAGTCCATTAAGGAGGCTTTATTGGCTGCT
 ACATGTTATTGCAAGAGACTTTGATGCCTTAAATGAACGCATCCAAAAGAGACAACAGAGCAGCGTGTCT
 CTTGAGGAACAAGAGAAAACAAGAAAGAGCTGAACGAGTGCGAAAATTACGAGAAGAAAGAAAACAAAATG
 AACAGGAGCAGGCTGAACTCGATGGAACCAAGTGGTCTGGCTGAGTTGGACCCAGAACCAACGAATCCTTT
 CCAGCCAATAGCATCTGTAATCATTGAGAATGAAGGAAAATTGAAAGAGAGAAAAAACCAAAAATG
 GAGAAAGACAGTGATGGCTGCCACCTGAAACATAAAATGGAGCATGAGCAATAGAGACACAAGGCCAGG
 TTAATCACAATGGCCAAAAAATAATGAATTTGGACTAGTAGAAAATTATAAGGAGGCATTAACACAGCA
 GTTAAAGAATGAAGATGAGACAGACCGGCCATCATTGGAATCAGCTAATGGTAAAAAGAAAATAAGAAA
 CTAAGAATGAAAAGGAACCACCGGTAGAACCCTTAATATAGATGACTGTGCTCCTGAGAGTCCAACGC
 CACCCCAACCCCTCCTCCTGTTGGCTGGGGAACCCCTAAAGTCACTAGACTTCCAAAACCTTGAGCCTCT
 TGGTGAAACACATCATAATGATTTCTATAGGAAGCCACTGCCTCCCCTGGCTGTGCCACAGCGACCTAAC
 AGTGATGCTCATGATGTGATCTCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC230067 representing NM_001174150
 Red=Cloning site Green=Tags(s)

MFSLMASCCGWFKRWREPVRKVTLMLVGLDNAGKTATAKGIQGEYPEDVAPTGVGFSKINLRQGFVETIF
 DLGGGIRIRGIWKNYYAESYGVIFVVDSSDEERMEETKEAMSEMLRHPRISGKPIILVLANKQDKEGALGE
 ADVIECLSLEKLVNEHKCLCQIEPCSAISGYGKKIDKSIKKGLYLLHVIARDFDALNERIQKETTEQRA
 LEEQEKQERAERVRKLRERKQNEQEAELDGTSGLAELDPEPTNPFQPIASVIEENEGKLEREKKNQKM
 EKDSGDGCHLKHKMEHEQIETQGQVNHNGQKNNEFGLVENYKEAL TQQLKNEDETDRPSLESANGKKTKK
 LRMKRNHRVEPLNIDDCAPESPTPPPPPVWGTPKVTRLPKLEPLGETHHNDFYRKPLPPLAVPQRPN
 SDAHDVIS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk8059_b11.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001174150

ORF Size: 1284 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001174150.2](#)

RefSeq ORF: 1287 bp

Locus ID: 200894


UniProt ID: [Q3SXY8](#)

Cytogenetics: 3q11.1-q11.2

MW: 49.1 kDa

Gene Summary:

This gene encodes a member of the ADP-ribosylation factor-like family. The encoded protein is a small GTPase that contains both N-terminal and C-terminal guanine nucleotide-binding motifs. This protein is localized in the cilia and plays a role in cilia formation and in maintenance of cilia. Mutations in this gene are the cause of Joubert syndrome 8. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]

Product images:


Circular map for RC230067