

Product datasheet for **RC230091L4V**

SRPX (NM_001170750) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SRPX
Synonyms:	DRS; ETX1; HEL-S-83p; SRPX1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001170750
ORF Size:	1332 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC230091).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001170750.1](#), [NP_001164221.1](#)

RefSeq ORF: 1335 bp

Locus ID: 8406

UniProt ID: [P78539](#)

Cytogenetics: Xp11.4

MW: 49.8 kDa



Gene Summary:

May be involved in phagocytosis during disk shedding, cell adhesion to cells other than the pigment epithelium or signal transduction.[UniProtKB/Swiss-Prot Function]