

Product datasheet for **RC230100L4V**

NIPAL4 (NM_001172292) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	NIPAL4
Synonyms:	ARCI6; ICHTHYIN; ICHYN; SLC57A6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001172292
ORF Size:	1341 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC230100).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001172292.1](#), [NP_001165763.1](#)

RefSeq ORF: 1344 bp

Locus ID: 348938

UniProt ID: [Q0D2K0](#)

Cytogenetics: 5q33.3

Protein Families: Transmembrane



MW: 48.4 kDa

Gene Summary: This gene likely encodes a membrane receptor. Mutations in this gene have been associated with autosomal recessive congenital ichthyosis. [provided by RefSeq, Feb 2010]