

## Product datasheet for **RC230122L4V**

### **GDPD2 (NM\_001171191) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	GDPD2 (NM_001171191) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GDPD2
Synonyms:	GDE3; OBDPF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001171191
ORF Size:	1380 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230122).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001171191.1</a> , <a href="#">NP_001164662.1</a>
RefSeq Size:	2072 bp
RefSeq ORF:	1383 bp
Locus ID:	54857
UniProt ID:	<a href="#">Q9HCC8</a>
Cytogenetics:	Xq13.1
Protein Families:	Transmembrane
MW:	52.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the glycerophosphodiester phosphodiesterase enzyme family. The encoded protein hydrolyzes glycerophosphoinositol to produce inositol 1-phosphate and glycerol. This protein may have a role in osteoblast differentiation and growth. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jan 2010]