

Product datasheet for **RC230140L4V**

GPR137 (NM_001170726) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GPR137 (NM_001170726) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR137
Synonyms:	C11orf4; GPR137A; TM7SF1L1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001170726
ORF Size:	1425 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230140).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001170726.1 , NP_001164197.1
RefSeq ORF:	1428 bp
Locus ID:	56834
UniProt ID:	Q96N19
Cytogenetics:	11q13.1
Protein Families:	Druggable Genome, Transmembrane
MW:	52.7 kDa



[View online »](#)

Gene Summary:

Lysosomal integral membrane protein that may regulate MTORC1 complex translocation to lysosomes (PubMed:31036939). May play a role in autophagy (PubMed:31036939). [UniProtKB/Swiss-Prot Function]