

Product datasheet for RC230152L3

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CTSA (NM_001167594) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: CTSA (NM_001167594) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: CTSA

Synonyms: GLB2; GSL; NGBE; PPCA; PPGB

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC230152).

Sgfl-Mlul

Sequence:

Restriction Sites: Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001167594

ORF Size: 1443 bp





CTSA (NM_001167594) Human Tagged Lenti ORF Clone - RC230152L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 001167594.1, NP 001161066.1

 RefSeq ORF:
 1446 bp

 Locus ID:
 5476

 UniProt ID:
 P10619

Cytogenetics:

Protein Families: Druggable Genome, Protease

20q13.12

Protein Pathways: Renin-angiotensin system

MW: 54.7 kDa

Gene Summary: This gene encodes a member of the peptidase S10 family of serine carboxypeptidases.


Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed to generate two chains that comprise the

heterodimeric active enzyme. This enzyme possesses deamidase, esterase and

carboxypeptidase activities and acts as a scaffold in the lysosomal multienzyme complex. Mutations in this gene are associated with galactosialidosis. [provided by RefSeq, Nov 2015]



Product images:



Circular map for RC230152L3