

## Product datasheet for **RC230160L3V**

### NEU4 (NM\_001167600) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NEU4 (NM_001167600) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NEU4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001167600
ORF Size:	1452 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230160).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001167600.1</a> , <a href="#">NP_001161072.1</a>
RefSeq ORF:	1455 bp
Locus ID:	129807
UniProt ID:	<a href="#">Q8WWR8</a>
Cytogenetics:	2q37.3
Protein Pathways:	Other glycan degradation, Sphingolipid metabolism
MW:	52 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to a family of glycohydrolytic enzymes, which remove terminal sialic acid residues from various sialo derivatives, such as glycoproteins, glycolipids, oligosaccharides, and gangliosides. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Nov 2009]