

## Product datasheet for **RC230177L4V**

### NEU4 (NM\_001167599) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NEU4 (NM_001167599) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NEU4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001167599
ORF Size:	1491 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230177).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001167599.2</a>
RefSeq Size:	2376 bp
RefSeq ORF:	1494 bp
Locus ID:	129807
UniProt ID:	<a href="#">Q8WWR8</a>
Cytogenetics:	2q37.3
Protein Pathways:	Other glycan degradation, Sphingolipid metabolism
MW:	53.1 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to a family of glycohydrolytic enzymes, which remove terminal sialic acid residues from various sialo derivatives, such as glycoproteins, glycolipids, oligosaccharides, and gangliosides. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Nov 2009]