

Product datasheet for RC230193L3

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

HCK (NM_001172129) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: HCK (NM_001172129) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: HCK

Synonyms: JTK9; p59Hck; p61Hck

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Sgfl-Mlul

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC230193).

Sequence:


<u>.</u>

Restriction Sites: Cloning Scheme:

Cloning sites used for ORF Shuttling:

Sgf 1 ORF Mlu I

GCG ATC GCC ATG ---//--- NNN ACG CGT ---



^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001172129

ORF Size: 1515 bp





HCK (NM_001172129) Human Tagged Lenti ORF Clone - RC230193L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 001172129.1

 RefSeq ORF:
 1518 bp

 Locus ID:
 3055

 UniProt ID:
 P08631

 Cytogenetics:
 20q11.21

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: Chemokine signaling pathway, Fc gamma R-mediated phagocytosis

MW: 57.8 kDa


Gene Summary: The protein encoded by this gene is a member of the Src family of tyrosine kinases. This

protein is primarily hemopoietic, particularly in cells of the myeloid and B-lymphoid lineages. It may help couple the Fc receptor to the activation of the respiratory burst. In addition, it may play a role in neutrophil migration and in the degranulation of neutrophils. Multiple isoforms with different subcellular distributions are produced due to both alternative splicing and the use of alternative translation initiation codons, including a non-AUG (CUG) codon. [provided


by RefSeq, Feb 2010]



Product images:



Circular map for RC230193L3



Double digestion of RC230193L3 using Sgfl and Mlul $\,$