

Product datasheet for **RC230215L4V**

TROVE2 (RO60) (NM_001173525) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TROVE2 (RO60) (NM_001173525) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RO60
Synonyms:	RORNPN; SSA2; TROVE2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001173525
ORF Size:	1554 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230215).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001173525.1 , NP_001166996.1
RefSeq ORF:	1557 bp
Locus ID:	6738
UniProt ID:	P10155
Cytogenetics:	1q31.2
Protein Families:	Transcription Factors
Protein Pathways:	Systemic lupus erythematosus
MW:	58.9 kDa



[View online »](#)

Gene Summary:

RNA-binding protein that binds to misfolded non-coding RNAs, pre-5S rRNA, and several small cytoplasmic RNA molecules known as Y RNAs. May stabilize some of these RNAs and protect them from degradation (PubMed:18056422). Binds to endogenous Alu retroelements which are induced by type I interferon and stimulate proinflammatory cytokine secretion. Regulates the expression of Alu retroelements as well as inflammatory genes (PubMed:26382853).[UniProtKB/Swiss-Prot Function]