

Product datasheet for **RC230217**

SIGLEC10 (NM_001171160) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SIGLEC10 (NM_001171160) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SIGLEC10
Synonyms:	PRO940; SIGLEC-10; SLG2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC230217 representing NM_001171160
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTACTGCCACTGCTGCTGCTCGCTGCTGGCGGGTCCCAGGCTATGGATGGGAGATTCTGGATAC
 GAGTGCAGGAGTCAAGTGGTGGCGGAGGGCCTGTGCATCTCTGTGCCCTGCTCTTTCTCTACCCCG
 ACAAGACTGGACAGGGTCTACCCAGCTTATGGCTACTGGTCAAAGCAGTGACTGAGACAACCAAGGGT
 GCTCCTGTGGCCACAACCACCAGAGTCGAGAGGTGAAATGAGCACCCGGGGCCGATTCCAGCTCACTG
 GGGATCCCGCCAAGGGAACTGCTCCTTGGTGATCAGAGACGCGCAGATGCAGGATGAGTACAGTACTT
 CTTTCGGTGGAGAGAGGAAGCTATGTGAGATAAATTCATGAACGATGGGTTCTTTCTAAAAGTAACA
 GCCCTGACTCAGAAGACCGTCCGACTCCGTGTGGCCTATGCCCCAGAGACCTTGTATCAGCATTTAC
 GTGACAACACGCCAGCCCTGGAGCCCAGCCCAGGGAAATGTCCCATACCTGGAAGCCAAAAGGCCA
 GTTCTGCGGCTCCTCTGTGCTGCTGACAGCCAGCCCCCTGCCACTGAGCTGGGTCTGCAGAACAGA
 GTCTCTCCTCGTCCATCCCTGGGGCCCTAGACCCCTGGGGCTGGAGCTGCCCGGGTGAAGGCTGGGG
 ATTCAGGGCGCTACACCTGCCGAGCGGAGAACAGGCTTGCTCCCAGCAGCGAGCCCTGGACCTCTCTGT
 GCAGTATCCTCCAGAACTGAGAGTGTGGTTCCCAAGCAAACAGGACAGTCTGGAAAACCTTGGG
 AACGGCACGTCTCTCCAGTACTGGAGGGCCAAAGCCTGTGCCTGGTCTGTGTCACACACAGCAGCCCC
 CAGCCAGGCTGAGCTGGACCCAGAGGGGACAGGTTCTGAGCCCCCTCCAGCCCTCAGACCCCGGGTCT
 GGAGCTGCCTCGGGTTCAAGTGGAGCACGAAGGAGGTTACCTGCCACGCTCGGCACCCACTGGGCTCC
 CAGCACGTCTCTCAGCCTCTCCGTGCACTAAGAAGGGACTCATCTAACGGCATTCTCAAACGGGAG
 CGTTTCTGGGAATCGGCATCAGGCTCTCTTTTCTCTGCTGGCCCTGATCATGAAGATTCTACC
 GAAGAGACGGACTCAGACAGAAACCCGAGGCCAGGTTCTCCCGCACAGCACGATCCTGGATTACATC
 AATGTGGTCCCGACGGCTGGCCCCCTGGCTCAGAAGCGGAATCAGAAAGCCACACCAAACAGTCTCGGA
 CCCCTCTCCACCAGGTCTCCCTCCCCAGAATCAAAGAAGAACCAGAAAAGCAGTATCAGTTGCCAG
 TTTCCAGAACCCAAATCATCCACTCAAGCCCCAGAATCCAGGAGAGCCAAGAGGAGCTCCATTATGCC
 ACGCTCAACTTCCAGGCGTCAGACCCAGGCTGAGGCCCGGATGCCAAGGGCACCCAGGCGGATTATG
 CAGAAGTCAAGTCCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC230217 representing NM_001171160
 Red=Cloning site Green=Tags(s)

MLLPLLLSSLLGGSQAMDGRFWIRVQESVMVPEGLCISVPCSFYPRQDWTGSPAYGYWFKAVTETTKG
 APVATNHQSREVMSTRGRFQLTGDPKGNCSLVIRDAQMQDESQYFRVERGSYVRYNFMNDGFFLKVT
 ALTQKTVRLRVAYAPRDLVISISRDNTPALEPQPQGNVPYLEAQKQFLRLLCAADSQPPATLSWVLQNR
 VLSSHPWGRPLGLELPGVKAGDSGRYTCRAENRLGSQQRALDLVQYPPENLRVMVSQANRTVLENLG
 NGTSLPVLQSLCLVCVTHSSPPARLSWTQRQVQLSPSQSPDGVLELPRVQVEHEGEFTCHARHPLGS
 QHVSLSLSVHYKGLISTAFSNGAFLGIGITALLFLCLALIIIMKILPKRRTQETETPRPRFSRHSTILDYI
 NVVPTAGPLAQKRQKATPNSPRTPPLPGAPSPESKKNQKQYQLPSFPEPKSSTQAPESQESQEELHYA
 TLNFPGVRPRPEARMPKGTQADYAEVKFQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk8065_b01.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_001171160

ORF Size: 1557 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001171160.1](#), [NP_001164631.1](#)

RefSeq ORF: 1559 bp

Locus ID: 89790


Cytogenetics: 19q13.41

Protein Families: Druggable Genome, Transmembrane

MW: 57.9 kDa

Gene Summary:

SIGLECs are members of the immunoglobulin superfamily that are expressed on the cell surface. Most SIGLECs have 1 or more cytoplasmic immune receptor tyrosine-based inhibitory motifs, or ITIMs. SIGLECs are typically expressed on cells of the innate immune system, with the exception of the B-cell expressed SIGLEC6 (MIM 604405).[supplied by OMIM, Jul 2002]

Product images:


Circular map for RC230217