

Product datasheet for **RC230239**

SHISA6 (NM_001173462) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	SHISA6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC230239 representing NM_001173462
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCGCTGCGGCGCTCCTGCTGCTGCTGCTGCTCTCGTGGAGTCCCTGGACCTGCTGCCAGCGTCC
 ACGGAGCCCGCGGCCGCGCCCAACCGGACCTGAGTGCAGGCGGCGCTGCCGTGCGGGGCGGAGGGC
 CGGGGGCGCCCTGGCACGGGGCGGCCGAGCTGAACGGCACCGCCCGGCGCCCGGAATCCCGGAGGCG
 GGAAGCCGGCGGGGCGAGCCCGCGGCGGCTGTGGCGGCGGCGGCGAGCGCGCCGTACCTACGAGACGT
 GCTGGGGCTACTACGACGTGAGCGGCCAGTACGACAAGGAGTTCGAGTGTAAACAGCGAGAGCGGCTA
 CCTGTACTGCTGCGGTACCTGCTACTACCGCTTCTGCTGCAAGAAGCGCCACGAGAAGCTGGACCAAGCGC
 CAGTGCACCAACTACCAGAGCCCGGTGTGGGTACAGACGCCAGCACCAAGGTGGTGTGCGCGGGGCCCCG
 AGAACAAGTACGACCCGGAAGGACAAGACCAACTTACCGTCTACATCACCTGCGGGGTATCGCCTT
 CGTCATCGTGGCGGCGTCTTCCCAAGGTCTCTACGACAAGGCCACCGCCCTCCACGGAGATGAAC
 ATCCACAGGGCTCTGGCTGACATCTTAAGACAACAGGGACCAATCCCCATAGCACACTGTGAAAGAGAAA
 CTATCTCGGCTATCGATACCTCTCCAAAGAGAACACGCCGGTCTAGATCGTCTCCAAAAACCACTACAC
 TCCTGTGCGTACGGCCAAGCAGACTCCAGGGCATTATGGGAAGGATGCTTACCGAAGTGGAGGACCTGAT
 TCCATAACTTCATCTCATCTGGATTGTACATTAGGAAGAGGACACCAAGGAGAAGCCACGGATGA
 ACAACATCCTGACATCAGCCACCGAGCCCTATGACCTCTCCTTCTCCCGCTCGTTCAGAACTTGGCCCA
 CCTGCCCCATCCTACGAGTCTGCAGTAAAGACCAACCCAGCAAGTATTCATCTCTGAAGAGGTAACC
 GACAAGGAGGTGACGAGTATTACATGAGACGGCGGCACCTGCCCGACCTGGCTGCCCGCGGCACCTCC
 CCCTCAATGTCATCCAGATGTCCCAACAGAAGCCGTTGCCAAGGGAACGACCCCGCGGCCATCCGGGC
 CATGTCCCAGGACAGGGTCTGTCCCGGATCGGGGCTGCCAGATGAGTTCAGCATGCCCTACGACCGC
 ATCCTGTCCGACGAGCAGCTGCTCTCCACGGAGCGCCTGCACTCCAGGACCCGCTGCTGTCCCGGAGC
 GGACGGCCTTTCCCGAGCAGTCTGTCCAGGGCCATCTCGCACACGGACGCTTTGTGTCCACACCCGT
 GCTGGACCGCTACCGCATGAGCAAGATGCACTCTCATCCAGTGCTCCAATAACTCATACGCCACCCTG
 GGCCAGAGCCAGACGGCAGCCAAGCGTCAATGCTTGCCTCAGCAGACACAACCGGTGGAGCAGCTGC
 ACTACATCCCGGGCCACACACCTGCTACACAGCCAGCAAGACCGAAGTACCCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTAA

Protein Sequence: >RC230239 representing NM_001173462
 Red=Cloning site Green=Tags(s)

MALRRLLLLLLSLESLLPSVHGARGRAANRTL SAGGA VGGRRAGGALARGGRELNGTARAPGIPEA
 GSRRGQPA AAVAAA AAVTYETCWGYDVSGQYDKEFE CNSESGLYCCGTCY YRFCKKRHEKLDQR
 QCTNYQSPVWVQTPSTKV VSPGPNKYDPEKDTNFTVYITCGVIAFVIVAGVFAKVS YDKAHRPPREMN
 IHRALADILRQQPIPIAH CERETISAIDTSPKENTPVRSSSKNHYTPVRTAKQTPGHYKDAYRSGGPD
 LHNFISSGFVTLGRGHTKEKPRMNNILTSATEPYDLFSR SFQNL AHLPPSYESA VKTNPSKYSSLKRLT
 DKEADEYYMRRRHLPDLAARGTLPLNVIQMSQQKPLPRERPRRPIRAMSQDRVLS PDRGLPDEF SMPYDR
 ILSDEQLLSTERLHSQDPLLSPERTAFPEQSLSRAISHTDVFVSTPVLDRYRMSKMHSHP SASNNSYATL
 GQSQTA AKRHAFASRRHNTVEQLHYIPGHHTCYTASKTEVT

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8065_g08.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001173462

ORF Size: 1596 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001173462.2](#)

RefSeq ORF: 1599 bp

Locus ID: 388336


UniProt ID: [Q6ZSJ9](#)

Cytogenetics: 17p12

MW: 59.6 kDa

Gene Summary: Involved in maintenance of high-frequency synaptic transmission at hippocampal CA3-CA1 synapses. Regulates AMPA-type glutamate receptor (AMPA) immobilization at postsynaptic density keeping the channels in an activated state in the presence of glutamate and preventing synaptic depression. May play a role in self-renewal and differentiation of spermatogonial stem cells by inhibiting canonical Wnt signaling pathway.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC230239