

Product datasheet for RC230241L3

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

STK23 (SRPK3) (NM_001170761) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: STK23 (SRPK3) (NM 001170761) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: STK23

Synonyms: MSSK-1; MSSK1; STK23

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)


E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC230241).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001170761

ORF Size: 1599 bp





STK23 (SRPK3) (NM_001170761) Human Tagged Lenti ORF Clone - RC230241L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 001170761.1, NP 001164232.1

RefSeq ORF: 1602 bp
Locus ID: 26576
UniProt ID: Q9UPE1
Cytogenetics: Xq28

Protein Families: Druggable Genome, Protein Kinase

MW: 59.4 kDa


Gene Summary: This gene encodes a protein kinase similar to a protein kinase which is specific for the SR

(serine/arginine-rich domain) family of splicing factors. A highly similar protein has been shown to play a role in muscle development in mice. Multiple transcript variants encoding

different isoforms have been found for this gene. [provided by RefSeq, Dec 2009]



Product images:



Circular map for RC230241L3