

## Product datasheet for **RC230250L3V**

### Asparagine synthetase (ASNS) (NM\_001178075) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Asparagine synthetase (ASNS) (NM_001178075) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Asparagine synthetase
Synonyms:	ASNSD; TS11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001178075
ORF Size:	1686 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230250).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001178075.1</a> , <a href="#">NP_001171546.1</a>
RefSeq Size:	2010 bp
RefSeq ORF:	1623 bp
Locus ID:	440
UniProt ID:	<a href="#">P08243</a>
Cytogenetics:	7q21.3
Protein Families:	Druggable Genome
Protein Pathways:	Alanine, aspartate and glutamate metabolism, Metabolic pathways, Nitrogen metabolism



[View online »](#)

**MW:** 64.4 kDa

**Gene Summary:** The protein encoded by this gene is involved in the synthesis of asparagine. This gene complements a mutation in the temperature-sensitive hamster mutant ts11, which blocks progression through the G1 phase of the cell cycle at nonpermissive temperature. Alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, May 2010]