

## Product datasheet for **RC230326L4V**

### STIM2 (NM\_001169117) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	STIM2 (NM_001169117) Human Tagged ORF Clone Lentiviral Particle
Symbol:	STIM2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001169117
ORF Size:	1797 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230326).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001169117.1</a> , <a href="#">NP_001162588.1</a>
RefSeq ORF:	1800 bp
Locus ID:	57620
UniProt ID:	<a href="#">Q9P246</a>
Cytogenetics:	4p15.2
MW:	68.5 kDa



[View online »](#)

**Gene Summary:**

This gene is a member of the stromal interaction molecule (STIM) family and likely arose, along with related family member STIM1, from a common ancestral gene. The encoded protein functions to regulate calcium concentrations in the cytosol and endoplasmic reticulum, and is involved in the activation of plasma membrane Orai Ca(2+) entry channels. This gene initiates translation from a non-AUG (UUG) start site. A signal peptide is cleaved from the resulting protein. Multiple transcript variants result from alternative splicing. [provided by RefSeq, Dec 2009]