

Product datasheet for **RC230332L4V**

IGSF10 (NM_001178146) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IGSF10 (NM_001178146) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IGSF10
Synonyms:	CMF608
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001178146
ORF Size:	1806 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230332).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001178146.1
RefSeq ORF:	1809 bp
Locus ID:	285313
UniProt ID:	Q6WRI0
Cytogenetics:	3q25.1
Protein Families:	Druggable Genome, Transmembrane
MW:	67.2 kDa


[View online »](#)

Gene Summary:

Involved in the control of early migration of neurons expressing gonadotropin-releasing hormone (GNRH neurons) (By similarity). May be involved in the maintenance of osteochondroprogenitor cells pool (By similarity).[UniProtKB/Swiss-Prot Function]