

Product datasheet for **RC230337**

SLC26A11 (NM_001166349) Human Tagged ORF Clone

Product data:

| | |
|--------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | SLC26A11 (NM_001166349) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | SLC26A11 |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |



[View online »](#)

ORF Nucleotide Sequence:

>RC230337 representing NM_001166349
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCTTCTTCGGTGACGGCGCTGGGTACAGCCAGGTCCTCTGGCCCCGGGATGGCCCCAGCGCCTGCT
 GCTGCTCCCTGCGGCCCTGCAGAGGAGGCTGCCCATCCTGGCGTGGCTGCCCAGCTACTCCTGCAGTG
 GCTGAAGATGGATTTCTGTCGCCGCTCTCAGTTGGCCTCACTGCCATCCCCAGGCGCTGGCCTATGCT
 GAAGTGGCTGGACTCCCAGCCAGTATGGCCTCTACTCTGCCTCATGGGCTGCTTCGTGTATTTCTTCC
 TGGGCACCTCCCGGATGTGACTCTGGGCCACCGCCATTATGTCCCTCCTGGTCTCCTTACACCTT
 CCATGAGCCCGCTACGCTGTGCTGGCCTTCTGTCCGGCTGCATCCAGCTGGCCATGGGGTCTGT
 CGTTTGGGTTCTGCTGGACTTCATTTCTACCCCGTCATTAAGGCTTCACTCTGCTGCTGCCGTCA
 CCATCGGCTTTGGACAGATCAAGAACCTGCTGGGACTACAGAACATCCCCAGGCCGTTCTTCTGCAGGT
 GTACCACACCTTCTCAGGATTGCAGAGACCAGGTTAGGTGACGCCGCTCCTGGGGCTGGTCTGCATGCTG
 CTGCTGCTGGTGTGAAGCTGATGCGGGACCACGTGCCTCCCGTCCACCCGAGATGCCCCCTGGTGTG
 GGCTCAGCCGTGGGCTGGTCTGGGCTGCCACGACAGCTCGCAACGCCCTGGTGGTCTCCTTCGACGCCCT
 GGTGCGTACTCCTTCGAGGTGACTGGATACCAGCCTTTCATCCTAACAGGGGAGACAGCTGAGGGGCTC
 CCTCCAGTCCGGATCCCGCCCTTCTCAGTGACCACAGCCAACGGGACGATCTCCTTACCAGGATGGTGC
 AGGACATGGGAGCCGGGCTGGCCGTGGTCCCTGATGGGCTCCTGGAGAGCATTGCGGTGGCCAAAGC
 CTTGCGATCTCAGAATAATTACCGCATCGATGCCAACCAGGAGCTGCTGGCCATCGGTCTACCAACATG
 TTGGGCTCCCTCGTCTCCTCTACCCGGTACAGGCAGCTTTGGACGGACAGCCGTGAACGCTCAGTCGG
 GGTGTGCACCCCGGGGGGCTGGTGACGGGAGTGTGGTGTGCTGTCTCTGGACTACCTGACCTC
 ACTGTTCTACTACATCCCAAGTCTGCCCTGGCTGCCGTATCATCATGGCCGTGGCCCGCTGTTTCGAC
 ACCAAGATCTTACAGGACGCTCTGGCGTGTAAAGAGGCTGGACCTGCTGCCCTGTGCGTACCTTCTG
 TGTGCTTCTGGGAGTGCAGTACGGCATCCTGGCCGGGCCCTGGTGTCTCTGCTCATGCTCCTGCACTC
 TGCAGCCAGGCTGAGACCAAGGTGTCAGAGGGGCCGTTCTGGTCTGCAGCCGGCCAGCGGCTGTCC
 TTCCCTGCCATGGAGGCTCTGCGGGAGGAGATCCTAAGCCGGGCCCTGGAAGTGTCCCGCCACGCTGCC
 TGGTCTGGAGTGCACCCATGTCTGCAGCATCGACTACACTGTGGTGTGGGACTCGCGAGCTCCTCCA
 GGACTTCCAGAAGCAGGGCGTGCCTGGCCTTTGTGGGCTGCAGGTCCCGTCTCCGTGCTGCTGCTG
 TCCGCTGACCTGAAGGGTTCCAGTACTTCTTACCCTGGAAGAAGCAGAGAAGCACCTGAGGCAGAAGC
 CAGGGACCCAGCCCTACAACATCAGAGAAGACTCCATTCTGGACCAAAGGTTGCCCTGCTCAAGGCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC230337 representing NM_001166349
 Red=Cloning site Green=Tags(s)

MPSSVTALGQARSSGPGMAPSACCCSPAALQRRPLILAWLPSYSLQWLKMDFVAGLSVGLTAIPQALAYA
 EVAGLPPQYGLYSAFMGCFVYFFLGTSRDVTLGPTAIMSLLVSFYTFHEPAYAVLLAFLSGCIQLAMGVL
 RLGFLLDIFISYPVIKGFSAAAVTIGFGQIKNLLGLQNIIPRPFLLQVYHTFLRIAETRVGDAVLGLVCM
 LLLVLKLMRDHVPVHPPEMPGVRLSRGLVWAATTARNALVVSFAALVAYSFEVTGYQPFILTGETAEG
 LPPVRIIPFVSVTANGTISFTEMVQDMGAGLAVVPLMGLLESIAVAKAFASQNNYRIDANQELLAIGLT
 NMLGSLVSSYPVTGSFGRYAVNAQSGVCTPAGGLVTGVLVLLSLDYLTSFYIIPKSALAAVIMAVAPLFD
 TKIFRTLWRVKRLDLLPLCVTFLLCFWEVQYIGILAGALVSLMLLHSAARPETKVSEGPVLLQPASGLS
 FPAMEALREEILSRALEVSPRCLVLECHTVCSIDYTVVLGLGELLQDFQKQGVVALAFVGLQVPLRVLL
 SADLKGQYFSTLEEAEKHLRQKPGTQPYNIREDSIDLQKVALLLKA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

| | |
|-------------------------------|---|
| ORF Size: | 1818 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001166349.1 , NP_001159821.1 |
| RefSeq ORF: | 1821 bp |
| Locus ID: | 284129 |
| UniProt ID: | Q86WA9 |
| Cytogenetics: | 17q25.3 |
| Protein Families: | Transmembrane |
| MW: | 65.7 kDa |
| Gene Summary: | This gene encodes a member of the solute linked carrier 26 family of anion exchangers. Members of this family of proteins are essential for numerous cellular functions including homeostasis and intracellular electrolyte balance. The encoded protein is a sodium independent sulfate transporter that is sensitive to the anion exchanger inhibitor 4,4'-diisothiocyanostilbene-2,2'-disulfonic acid. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Oct 2009] |