

Product datasheet for **RC230342**

BRIT1 (MCPH1) (NM_001172574) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BRIT1 (MCPH1) (NM_001172574) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	BRIT1
Synonyms:	BRIT1; MCT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC230342 representing NM_001172574
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**C

ATGGCGGCCCATCCTGAAAGATGTAGTGGCCTATGTTGAAGTGTGGTCATCCAATGGAACAGAAAATT
 ATTCAAAGACATTTACAACACAGCTTGTGGATATGGGGGCAAAGGTTTCAAAAACCTTTAAACAAACAAGT
 AACTCACGTTATCTTCAAAGATGGCTACCAGAGCACTTGGGACAAAGCTCAGAAGAGAGCGTAAAGCTC
 GTTTCGGTGCTCTGGGTTGAAAAATGCAGGACAGCTGGAGCACACATTGATGAATCATTGTTCCCTGCAG
 CTAATATGAATGAACACTTATCAAGCCTAATTAACAAAAACGTAATGTATGCAGCCCAAGGATTTTAA
 TTTTAAACACCAGAAAATGATAAGAGATTTTCAAGAAAATTTGAGAAAATGGCTAAAGAGCTACAAAGG
 CAAAAACAAATCTAGATGATGATGTACCTATTCTCTTATTTGAATCTAATGGTTCATTAATATATACTC
 CCACAATTGAAATTAATAGTAGTCACCACAGCGCAATGGAGAAGAGATTACAAGAGATGAAGGAGAAAAAG
 GGAAAATCTTTCCCCACCTCTTCCCAATGATTACGAGTCTCATGATAATCCAAGTAACCTCTGTGT
 GAAGCACCTTTGAACATTTACGTGATACCTTTGTGTTTCAATGAATACTTTGCTGGTGGCTTACACTCAT
 CTTTTGATGATCTTTGTGAAACTCAGGATGTGAAATCAGGAAAGGAAGTTGGAAGGATCCATTATGA
 CATTAAAGTGATGTGTGATTTTCTTCACTTGTATTGAAAGCAAATAATATTCAATCATCACCATCTTTC
 ACTCACCTCGATAAATCAAGTCCTCAGAAATTTCTGAGTAATCTTTCAAAGGAAGAAATAAACTTGCAAA
 GAAATATTGCAGGTAAAGTAGTCACCCCTCACCAAAAGCAGGCTGCAGGTATGTCTCAGGAGACGTTTGA
 AGAGAAGTATCGTTTGTCTCCTACCTTATCTTCAACAAAAGGCCACCTTTTGATACATTCAGACCCAGG
 AGTTCTCAGTAAAGAGAAAAAGAGTATCATGGCTCCCATTCACCTCCGAAGGAAAAATGCAAGAGAA
 AGAGGAGCACCAGGAGATCTATCATGCCGAGGCTGCAGCTGTGCAGGTCCGAAGGCAGGCTGCAGCACGT
 GCGGGGACCTGCCCTGGAGGCTCTTAGCTGTGGGAGTCTTCATATGATGACTATTTTTCACCTGATAAT
 CTTAAGGAAAGGTATTACAGAAATCTTCTCCTGAATCTCAGCTGCCATCAAGCCCTGCTCAGTTGAGCT
 GCAGAAGTCTTTCTAAGAAGGAGAGAAACAGCATATTTGAAATGTCTGATTTTCTGCTTGGCAAAAA
 AACCAGAACAGTTGACATTACCAATTTACAGCAAAAACCATCTCCAGTCTCGGAAAACCTGGAATGGT
 GAAGGCCGTGCAACTTCGAGTTGCGTGACTTCTGCCCTGAAGAAGCCCTAAGGTGTTGTAGACAGGCTG
 GGAAAGAAGACGCATGCCAGAGGGAATGGCTTTTCTTACACCATTGAGGACCTGCTCTTCCAAAAGG
 ACATGATGATGATTTAACTCCTTTGGAAGGAAGCCTTGAAGAAATGAAAGAAGCGGTTGGTCTGAAAAGC
 ACACAGAACAAAGGTACCACTTCCAAAATATCAAACTCCTCTGAAGGCGAAGCCAGAGTGAACATGAGC
 CATGTTTTATAGTTGACTGTAACATGGAGACGTCTACAGAAGAGAAGGAAAACCTACCCGGAGGATACAG
 TGAAGATATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC230342 representing NM_001172574
 Red=Cloning site Green=Tags(s)

MAAPILKDVVAYEVWSSNGTENYSKFTTQLVDMGAKVSKTFNKQVTHVIFKDGYSQSTWDKAQKRGVKL
 VSVLWVEKCRTAGAHIDESLFPANMNEHLSSLIKKKRCKMQPKDFNFKTPENDKRFQKKFEKMAKELQR
 QKTNLDDVPILLFESNGSLIYPTIEINSSHSAAMEKRLQEMKEKRENLSPTSSQMIQQSHDNPSNSLC
 EAPLNISRDTLCSDEYFAGGLHSSFDDLCGNSGCGNQRKLEGSINDIKSDVCISLVLKANNIHSSPSF
 THLDKSSPQKFLSNLSKEEINLQRNIAGKVTPHQKQAAGMSQETFEKYRLSPTLSSTKGHLIIHSRPR
 SSSVKKRVSHGSHSPPKEKCKRKRSTRRSIMPRLLQLCRSEGRQLQHVAGPALEALSCGESSYDDYFSPDN
 LKERYSENLPPESQLPSSPAQLSCRSLSKKERTSIFEMSDFCVGGKTRTVDTNFTAKTISSPRKTGNG
 EGRATSSCVTSAPEEALRCCRQAGKEDACPEGNGFSYTIEDPALPKGHDDDLTPLEGSLEEMKEAVGLKS
 TQNKGTTSKISNSSEGEAQSEHEPCFIVDCNMETSTEEKENLPGGYS GSM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001172574

ORF Size: 1830 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001172574.1](#), [NP_001166045.1](#)

RefSeq ORF: 1833 bp

Locus ID: 79648


UniProt ID: Q8NEM0

Cytogenetics: 8p23.1

MW: 68.3 kDa

Gene Summary: This gene encodes a DNA damage response protein. The encoded protein may play a role in G2/M checkpoint arrest via maintenance of inhibitory phosphorylation of cyclin-dependent kinase 1. Mutations in this gene have been associated with primary autosomal recessive microcephaly 1 and premature chromosome condensation syndrome. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2010]

Product images:



Circular map for RC230342