

Product datasheet for **RC230377**

SYK (NM_001174167) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SYK (NM_001174167) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SYK
Synonyms:	IMD82; p72-Syk
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RC230377 representing NM_001174167
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCCAGCAGCGGCATGGCTGACAGCGCCAACCACTGCCCTTCTTTTTCGGCAACATCACCCGGGAGG
AGGCAGAAGATTACCTGGTCCAGGGGGCATGAGTGATGGGCTTTATTTGCTGCGCCAGAGCCGCAACTA
CCTGGGTGGCTTCGCCCTGTCCGTGGCCACGGGAGGAAGGCACACCACTACACCATCGAGCGGGAGCTG
AATGGCACCTACGCCATCGCCGGTGGCAGGACCCATGCCAGCCCCGCGACCTCTGCCACTACCACTCCC
AGGAGTCTGATGGCCTGGTCTGCCTCCTCAAGAAGCCCTCAACCGGCCCAAGGGTGCAGCCCAAGAC
TGGGCCCTTTGAGGATTTGAAGGAAAACCTCATCAGGGAATATGTGAAGCAGACATGGAACCTGCAGGGT
CAGGCTCTGGAGCAGGCCATCATCAGTCAGAAGCCTCAGCTGGAGAAGCTGATCGCTACCACAGCCCATG
AAAAATGCCTTGGTTCCATGAAAAATCTCTCGGGAAGAATCTGAGCAAATTCCTGATAGGATCAA
GACAAATGAAAGTTCCATCGAGCCAGAGACAACAACGGCTCTACGCCCTGTGCCTGCTGCACGAA
GGGAAGGTGCTGCACTATCGCATCGACAAGACAGGGAAGCTCTCCATCCCCGAGGAAAGAAGT
TCGACACGCTCTGGCAGCTAGTCGAGCATTATTCTTATAAAGCAGATGGTTTGTTAAGAGTCTTACTGT
CCCATGTCAAAAAATCGGCACACAGGAAATGTTAATTTTGGAGGCCCTCCACAACCTCCAGGTTCCCAT
CCTGCGACTTGGTCAGCGGGTGAATAATCTCAAGAATCAAATCATACTCCTTCCCAAAGCTGGCCACA
GAAAGTCTCCCTGCCAAGGGAACCGGCAAGAGAGTACTGTGTCATTCAATCCGTATGAGCCAGA
TGCACCCCTGGGCTGCAGACAAAGGCCCCAGAGAGAAGCCCTACCATGGACACAGAGGTGTACGAGAGC
CCCTACGCGGACCCTGAGGAGATCAGGCCAAGGAGGTTTACCTGGACCGAAAGCTGCTGACGCTGGAAG
ACAAAGA
AAGTCTGGCTCTGGTAATTTTGGAACTGTGAAAAAGGGCTACTACCAAATGAAAAAAGTTGTGAA
AACCGTGGCTGTGAAAAACTGAAAAACGAGGCCAATGACCCCGCTCTTAAAGATGAGTTATTAGCAGAA
GCAAATGTCATGCAGCAGCTGGACAACCCGTACATCGTGCATGATCGGGATATGCGAGGCCGAGTCTCT
GGATGCTAGTTATGGAGATGGCAGA
AAGTGGTCCCTCAATAAGTATTTGCAGCAGAACAGACATGTCAA
GGATAAGAACATCATAGA
AAGTGGTTCATCAGGTTCCATGGGCATGAAGTACTTGGAGGAGAGCAATTTT
GTGCACAGAGATCTGGCTGCAAGAAATGTGTTGCTAGTTACCAACATTATGCCAAGATCAGTGATTTG
GACTCTCCAAAGCACTGCGTGTGATGAAA
AAGTACTACAAGGCCAGACCCATGGAAAGTGGCCTGTCAA
GTGGTACGCTCCGAATGCATCAACTACTACAAGTTCTCCAGCAAAGCGATGTCTGGAGCTTTGGAGTG
TTGATGTGGGAAGCATTCTCTATGGGCAGAAGCCATATCGAGGGATGAAAGGAAGTGAAGTCAACGCTA
TGTTAGAGAAAGGAGAGCGGATGGGGTGCCTGCAGGGTGTCCAAGAGAGATGTACGATCTCATGAATCT
GTGCTGGACATACGATGTGAAAAACAGGCCCGATTTCGCAGCAGTGGAACTGCGGCTGCGCAATTACTAC
TATGACGTGGTGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC230377 representing NM_001174167
 Red=Cloning site Green=Tags(s)


MASSGMADSANHLPPFFGNITREEAEDYLVQGGMSDGLYLLRQSRNYLGGFALSVAHGRKAHHTIEREL
 NGTYAIAGGRTHASPADLCHYHSQESDGLVCLLKKPFNRQGVQPKTGPFDLKENLIREYVKTWNLQG
 QALEQAIISQKPQLEKLIATTAHEKMPWFHGKISREESEQIVLIGSKTNGKFLIRARDNNGSYALCLLHE
 GKVLHYRIDKDKTGKLSIPEGKKFDLWQLVEHYSYKADGLLRVLTVPQKIGTQGNVNFVGGRRPQLPGSH
 PATWSAGGIISRIKSYSFPKPGHRKSSPAQGNRQESTVSFNPEPELAPWAADKGPQREALPMDTEVYES
 PYADPEEIRPKEYLDRKLLTLEDKELGSGNFGTVKKGYQMKVVKTVAVKILKNEANDPALKDELLAE
 ANVMQQLDNPIYVRMIGICEAESWMLVMEAEGLPNKYLQQRNHVKDKNIIELVHQVSMGMKYLEESNF
 VHRDLAARNVLLVTQHYAKISDFGLSKALRADENYYKAQTHGKWPVKWYAPECINYYKFSKSDVWSFGV
 LMWEAFSYGQKPYRGMKGSEVTAMLEKGERMGCPAGCPREMYDLMNLCTYDVENRPGFAAVELRLRNY
 YDVVN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

Plasmid Map:


ACCN: NM_001174167

ORF Size: 1905 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001174167.3](#)

RefSeq ORF:	1908 bp
Locus ID:	6850
UniProt ID:	P43405
Cytogenetics:	9q22.2
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	B cell receptor signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Natural killer cell mediated cytotoxicity
MW:	72.5 kDa
Gene Summary:	<p>This gene encodes a member of the family of non-receptor type Tyr protein kinases. This protein is widely expressed in hematopoietic cells and is involved in coupling activated immunoreceptors to downstream signaling events that mediate diverse cellular responses, including proliferation, differentiation, and phagocytosis. It is thought to be a modulator of epithelial cell growth and a potential tumour suppressor in human breast carcinomas. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2010]</p>