

Product datasheet for **RC230379**

PREPL (NM_001171617) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PREPL (NM_001171617) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PREPL
Synonyms:	CMS22
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC230379 representing NM_001171617
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGATGCATTTGAAAAAGTGAGAACAAAATTAGAAACACAGCCACAAGAAGAATATGAAATCATCAATG
TGGAAAGTTAAACATGGTGGTTTTGTTTATTACCAAGAAGTTGTTGCTTGGTTTCGTTCCAAAGATGAAGA
AGCAGACAATGATAATTATGAAGTTTTATTCAATTTGGAGGAACCTAAGTTAGACCAGCCCTTCATTGAT
TGTATCAGAGTTGCTCCAGATGAAAAATATGTGGCTGCCAAGATAAGAAGTGAAGATTCTGAAGCATCTA
CCTGTGTAATTATAAAGCTCAGCGATCAGCCCGTAATGGAAGCTTCTTTCCCGAATGTGTCCAGTTTTGA
ATGGGTAAGGACGAGGAAGATGAAGATGTTTTATTCTACACCTCCAGAGGAACCTTCGCTGTCATGAC
GTATATCGAGCCACTTTTGGTGATAACAAACGTAATGAACGCTTTTACACAGAAAAAGACCCAAGCTACT
TTGTTTTCTTTATCTTACAAAAGACAGTCGTTTCTCACCATAAATATTATGAACAAGACTACTTCTGA
AGTGTGGTTGATAGATGGCCTGAGCCCTGGGACCCACCAGTACTTATCCAGAAGCGAATACATGGGGTC
CTTTACTATGTTGAACACAGAGATGATGAATTATACATTCCTCACTAATGTTGGAGAACCTACAGAATTTA
AGCTAATGAGAACAGCGCTGATACCCCTGCAATTATGAATTGGGATTTATTTTTACAATGAAGAGAAA
TACAAAAGTGATAGACTTGGACATGTTAAGGATCACTGTGTTCTATTTCTGAAGCACAGCAATCTCCTT
TATGTTAATGTGATTGGTCTGGCTGATGATTCAGTTCCGGTCTCTAAAGCTCCCTCCTTGGGCTGTGGAT
TCATAATGGATACAAATCTGACCCAAAGAAGTCCCTTTCACTTTGCTCTCCAATACGTCCCCAAA
ATATTACACATACAAGTTTGCAGAAGGCAAACTGTTTGGGAACTGGGCATGAAGACCAATCACAAAG
AACTGACTCTGAGGACTTGCAGAAGAAACCTCTCTTGGTACATGTATATGGAGCTTATGGAATGGATTT
GAAAATGAATTTACGGCTGAGAGGCGGGTCTGGTGGATGATGGATGGATATTAGCATACTGCCATGTT
CGAGGTGGTGGTGAAGTAGGCCTCCAGTGGCACGCTGATGGCCGCTAACTAAAAAACTCAATGGCCTTG
CTGATTTAGAGGCTTGCATTAAGACGCTTCATGGCCAAGGCTTTTCTCAGCCAAGTCTAACAACCTGAC
TGCTTTCAGTGTGGAGGGTGTGCTGCAGGAGCATTGTGTAATTCTAATCCAGAGCTGGTGAAGCGGTG
ACTTTGGAGGCACCTTTCTTGGATGTTCTCAACACCATGATGGACTACACTTCTCTGACATTAGAAG
AATTAGAAGAATGGGGAATCCTTCATCTGATGAAAAACACAAGAAGTACATAAAAACGTTACTGTCCCTA
TCAAAATATTAACCTCAGCATTATCCTTCAATTCACATAACGGCATATGAAAACGATGAACGGTACCT
CTGAAAGGAATTGTAAGTTATACTGAGAACTCAAGGAAGCCATCGCGGAGCATGCTAAGGACACAGGTG
AAGGCTATCAGACCCCTAATATTATTCTAGATATTCAGCCTGGAGGCAATCATGTAATTGAGGATTCTCA
CAAAAAGATTACAGCCAAATTAATTCCTGTACGAGGAAGTGGACTTGACAGCACAGTGTTCGAG
GATCTTAAGAAATACCTGAAATTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA

Protein Sequence: >RC230379 representing NM_001171617
Red=Cloning site Green=Tags(s)

```


MDAFEKVRTKLETQPQEEYEIINVEVKHGGFVYYQEGCCLVRSKDEEADNDNYEVLFNLEELKLDQPFID
CIRVAPDEKYVAAKIRTEDSEASTCVIIKLSQPVMEASFNVSSFVWVWVDEEDEDVLFYTFQRNLRCHD
VYRATFGDNKRNERFYTEKDPSTYFVFLYLTKDSRFLTIMINMKTSEVWLIDGLSPWDDPPVLIQKRIHGV
LYYVEHRDDELYILTNVGEPTFEKLMRTAADTPAIMNWDLFFTMKRNTKVIDLDMFKDHCVLFLKHSNLL
YVNVIGLADDSVRSCLKPPWACGFIMDTNSDPKNCFQLCSPIRPPKYTYKFAEGKLFEEFGHEDPITK
TSRVLRLLEAKSKDGKLVPMTVFHKTDSEDLQKKPLLHVHYGAYGMDLKMNFRPERRVLVDDGWILAYCHV
RGGGELGLQWHADGRLTKKLNGLADLEACIKTLHGQGFSSQPSLTTLTAFSAGGLVLAGALCNSNPVLRVAV
TLEAPFLDVLNTMDDTTLPLTLEEEWGNPSSDEKHKNYIKRYCPYQNIKPQHYSIHITAYENDERVP
LKGIVSYTEKLEAIAEHAKDTGEGYQTPNIILDIQPGGNHVIEDSHKKITAQIKFLYEELGLDSTSVFE
DLKKYLKF
    
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8070_h01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001171617

ORF Size: 1914 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001171617.1](#), [NP_001165088.1](#)

RefSeq Size: 6064 bp

RefSeq ORF: 1917 bp

Locus ID: 9581

UniProt ID: [Q4J6C6](#)


Cytogenetics: 2p21

Protein Families: Druggable Genome, Protease

MW: 73.4 kDa

Gene Summary: The protein encoded by this gene belongs to the prolyl oligopeptidase subfamily of serine peptidases. Mutations in this gene have been associated with hypotonia-cystinuria syndrome, also known as the 2p21 deletion syndrome. Several alternatively spliced transcript variants encoding either the same or different isoforms have been described for this gene.[provided by RefSeq, Jan 2010]

Product images:



Circular map for RC230379