

## Product datasheet for **RC230384L2V**

### Zinc finger protein 668 (ZNF668) (NM\_001172669) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | Zinc finger protein 668 (ZNF668) (NM_001172669) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | Zinc finger protein 668  |
| Mammalian Cell Selection: | None   |
| Vector:                   | pLenti-C-mGFP (PS100071)   |
| Tag:                      | mGFP   |
| ACCN:                     | NM_001172669   |
| ORF Size:                 | 1926 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC230384).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001172669.1</a> , <a href="#">NP_001166140.1</a>  |
| RefSeq ORF:               | 1929 bp  |
| Locus ID:                 | 79759  |
| UniProt ID:               | <a href="#">Q96K58</a>   |
| Cytogenetics:             | 16p11.2  |
| MW:                       | 71 kDa   |
| Gene Summary:             | May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]  |



[View online »](#)