

## Product datasheet for **RC230406L4V**

### CPXM (CPXM1) (NM\_001184699) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CPXM (CPXM1) (NM_001184699) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CPXM
Synonyms:	CPX1; CPXM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001184699
ORF Size:	1980 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230406).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001184699.1</a>
RefSeq ORF:	1983 bp
Locus ID:	56265
UniProt ID:	<a href="#">Q96SM3</a>
Cytogenetics:	20p13
Protein Families:	Druggable Genome, Protease, Secreted Protein
MW:	73.9 kDa



[View online »](#)

**Gene Summary:**

This gene likely encodes a member of the carboxypeptidase family of proteins. Cloning of a comparable locus in mouse indicates that the encoded protein contains a discoidin domain and a carboxypeptidase domain, but the protein appears to lack residues necessary for carboxypeptidase activity.[provided by RefSeq, May 2010]