

Product datasheet for **RC230407**

ZNF185 (NM_001178109) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF185 (NM_001178109) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ZNF185
Synonyms:	SCELL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC230407 representing NM_001178109
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGTATCTCAGCTCTTGGAGGCCGACCAAAGGGAAGCCTCTGCCACCAGGCGAGGAGGAGCGCAATA
ACGTTCTCAAGCAGATGAAAGTGCGAACACGCTGAAGGGGGACAAGAGCTGGATTACCAAGCAGGATGA
ATCGGAGGGTCGCACCATAGAGCTGCCCTCAGGCCGAGTCGCGCCACATCCTTTTCATCAGCTGGGGAG
GTTCCGAAGCCTAGGCCTCCGAGACAAGGGCTCCCACTGGCTACATCATCCGGGGAGTGTTCACCAAGC
CCATAGACAGCTTCCAGCCCCAGCAGCAGTCCCAAGGCCAACGGGACTCCAAAAGTGCTGCCAG
TCTGGTGAGAACGCTAACGCTGGTCTCCCGCCCTCTCTCTGGCTACAAGATGACCACTGAGGAT
TACAAGAAGCTGGCACCTACAATATCAGGCGCAGCTCTACATCAGGGGACACCGAGGAGGAGGAGG
AGGAGGTGGTGCATTCTCTCAGATGAACAGAAACGGAGGTGAGGGCTGCAAGCGGTGTTCTGAGGAG
GACAGCTCCCGGGAGCACTCTACGTCTGTGAGCGGCAAGAAGAGCACTGGTCTACCCAGGAGACA
CAGGCACCGTTTATCGCGAAGAGGGTGGAGGTGGTGAAGAGGACGGGCCCTCTGAGAAGAGCCAGGACC
CACCTGCTCTGGCAAGATCCACTCTGGCTCAAACAGTGCAGGATGGAGGCAGGACCAAAGCGTCTCGGGC
AATTTGGATCGAGTGCCTGCCAAGTATGCCTAGCCCCGCTGGGAGCCAGGAGCTCAGGTCTTCCCAGGC
AACAAAGACAAGGAGGCCCTGCTCCAGAGAGCTCCAGAGGGACTTGGCTGGTGAAGGAGGCTTTCAGGG
CCCCAACACAGATGCTGCAAGGTCAAGTGCACAGTTGAGTGTGGCAATGTGGGATCCGGAGCCACGGG
CTCCCGCCTGAAGGCTTGGCTGCAGTAGACATCGGCTCCGAGAGAGGAAGCTCCAGTCCACTTTCAGTC
TCTGCTGTCCCTGCTGATAGGAAGAGCAACAGCACAGCAGCCAGGAGGATGCAAAGGCAGACCCAAAGG
GGCCTTGGCTGATTATGAGGGGAAGGATGTGGCCACCAGGGTCGGAGAGGCCTGGCAGGAGAGGCCTGG
AGTCCAAGAGGTGGCCAAGGAGACCCAGCTGTACCCGCTCAGCAACCTGCAGATCCCAGCACCCAGAG
CGGCAGAGCAGCCCCAGCGGATCTGAGCAACTTGTGACAGGAGAGTGTGGCAGCAGCGTGTGACTG
ATTTTGAGGGGAAGGATGTGGCCACCAAGGTCGGAGAGGCCTGGCAGGACAGGCCTGGAGCCCCAAGAGG
TGGCCAAGGAGACCCAGCTGTACCCACTCAGCAACCTGCAGATCCAGTACCCAGAACAGCAGAACAGC
CCCAGCGGATCTGAGCAATTCGTCAGACGAGAGAGCTGCACCAGCAGGGTGAAGGACCCCTCGAGCTGCA
TGGTCACTGTTACTGTCAGTCCACATCTGAGCAGCCTCACATTTATATCCAGCCCCGCAAGTGAATT
GGACTCCAGCTCTACCACAAAGGGATTCTCTTCGTGAAGGAGTACGTGAATGCTAGTGAAGTGTCTTCT
GGGAAGCCAGTATCTGCACGCTATAGCAACGTCAGCAGCATTGAGGACTCATTGCCATGGAGAAGAAGC
CTCCATGTGGCAGCACTCCATACTCTGAGAGGACAACCTGGAGGGATCTGTACTTACTGCAACCGTGAGAT
CCGAGACTGTCCAAAGATTACCTAGAACATCTTGGTATCTGCTGCCATGAATATTGCTTTAAGTGTGGG
ATTTGAGTAAACCGATGGGCGATCTCCTGGATCAGATCTTCAATCACCGTGACACCATTCACTGTGGGA
AATGCTATGAGAAGCTCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC230407 representing NM_001178109
 Red=Cloning site Green=Tags(s)


MSISALGGRTKGKPLPPGEEERNNVLKQMKVRTTLKGDKSWITKQDESEGRTEI LPSGRSRATSFSSAGE
 VPKPRPPSTRAPTYIIRGVFTKPIDSSSQPQQQFPKANGTPKSAASLVRTANAGPPRPSSSGYKMTTED
 YKKLAPYNI RRSSTSGDTEEEEEEEVVPFSSDEQKRRSEAA SGVLRRTPREHSYVLSA AKKSTGPTQET
 QAPFIAKRVEVVEEDGPSEKSQDPPALARSTPGSNSADGGRTKASRAIWIECLPSMPSPAGSQELRSPG
 NKDKEAPCSRELQRDLAGEEAFRAPNTDAARSSAQLSDGNVSGSATGSRPEGLAAVDIGSERGSSSATS
 SAVPADRKSNTAAQEDAKADPKGALADYEGKDVATRVEAWQERP GAPRGGQGDPAVPAQPADPSTPE
 RQSSPSGSEQLVRESCGSSVL TDFEGKDVATKVGEAWQDRPGAPRGGQGDPAVPTQQPADPSTPEQQNS
 PSGSEQFVRRESCTSRVSPSSCMVTVTATSEQPHIYIPAPASELDSSTTKGILFVKEYVNA SEVSS
 GKPVSARYSNVSSIEDSFAMEKKPPCGSTPY SERTTGGICTYCNREIRDCPKITLEHLGICCHEYCFKCG
 ICSKPMGDLLDQIFIHRDTIHCGKCYEKL F

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



ACCN: NM_001178109

ORF Size: 1980 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001178109.1](#), [NP_001171580.1](#)

RefSeq ORF: 1983 bp

Locus ID: 7739


UniProt ID: [O15231](#)

Cytogenetics: Xq28

MW: 70.8 kDa

Gene Summary: Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This gene encodes a LIM-domain zinc finger protein. The LIM domain is composed of two contiguous zinc finger domains, separated by a two-amino acid residue hydrophobic linker. The LIM domain mediates protein:protein interactions. Multiple alternatively spliced transcript variants encoding different isoforms have been identified.[provided by RefSeq, May 2010]

Product images:



Circular map for RC230407