

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for RC230416L3V

Kaiso (ZBTB33) (NM_001184742) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Kaiso (ZBTB33) (NM_001184742) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Kaiso
Synonyms:	ZNF-kaiso; ZNF348
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001184742
ORF Size:	2016 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230416).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 001184742.1</u>
RefSeq ORF:	2019 bp
Locus ID:	10009
UniProt ID:	<u>Q86T24</u>
Cytogenetics:	Xq24
MW:	74.9 kDa



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Kaiso (ZBTB33) (NM_001184742) Human Tagged ORF Clone Lentiviral Particle – RC230416L3V

Gene Summary:This gene encodes a transcriptional regulator with bimodal DNA-binding specificity, which
binds to methylated CGCG and also to the non-methylated consensus KAISO-binding site
TCCTGCNA. The protein contains an N-terminal POZ/BTB domain and 3 C-terminal zinc finger
motifs. It recruits the N-CoR repressor complex to promote histone deacetylation and the
formation of repressive chromatin structures in target gene promoters. It may contribute to
the repression of target genes of the Wnt signaling pathway, and may also activate
transcription of a subset of target genes by the recruitment of catenin delta-2 (CTNND2). Its
interaction with catenin delta-1 (CTNND1) inhibits binding to both methylated and non-
methylated DNA. It also interacts directly with the nuclear import receptor Importin-α2
(also known as karyopherin alpha2 or RAG cohort 1), which may mediate nuclear import of
this protein. Alternatively spliced transcript variants encoding the same protein have been
identified.[provided by RefSeq, May 2010]

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US