

## Product datasheet for **RC230419L3V**

### VIT (NM\_001177969) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	VIT (NM_001177969) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VIT
Synonyms:	VIT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001177969
ORF Size:	2034 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230419).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001177969.1</a> , <a href="#">NP_001171440.1</a>
RefSeq ORF:	2037 bp
Locus ID:	5212
UniProt ID:	<a href="#">Q6UXI7</a>
Cytogenetics:	2p22.2
Protein Families:	Secreted Protein
MW:	74.4 kDa



[View online »](#)

**Gene Summary:**

This gene encodes an extracellular matrix (ECM) protein. The protein may be associated with cell adhesion and migration. High levels of expression of the protein in specific parts of the brain suggest its likely role in neural development. [provided by RefSeq, Jun 2016]