

Product datasheet for **RC230431**

SLC34A2 (NM_001177999) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC34A2 (NM_001177999) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SLC34A2
Synonyms:	NAPI-3B; NAPI-IIb; NPTIIb; PULAM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC230431 representing NM_001177999
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCTCCCTGGCTGAATTGGGAGATGCCAGCCCAACCCCGATAAGTACCTCGAAGGGGCCGAGGTC
AGCAGCCCACTGCCCTGATAAAAGCAAAGAGACCAACAAAAAATAACACTGAGGCACCTGTAACCAAGAT
TGAACCTCTGCCGTCTACTCCACGGCTACACTGATAGATGAGCCCACTGAGGTGGATGACCCCTGGAAC
CTACCCACTCTTCAGGACTCGGGGATCAAGTGGTCAGAGAGAGACACCAAGGGAAGATTCTCTGTTTCT
TCCAAGGGATTGGGAGATTGATTTTACTTCTCGGATTTCTACTTTTTTCGTGTGCTCCCTGGATATTCT
TAGTAGCGCTTCCAGCTGGTTGGAGGAAAAATGGCAGGACAGTTCTTCAGCAACAGCTCTATTATGTCC
AACCTTTGTTGGGGCTGGTGATCGGGGTGCTGGTGACCGTCTTGGTGCAGAGCTCCAGCACCTCAACGT
CCATCGTTGTGAGCATGGTGTCTTTCATTGCTCACTGTTGCGGCTGCCATCCCATTATCATGGGGC
CAACATTGGAACGTCAATCACCAACACTATTGTTGCGCTCATGCAGGTGGGAGATCGGAGTGAGTTCAGA
AGAGCTTTTGCAGGAGCCACTGTCCATGACTTCTCAACTGGCTGTCCGTGTTGGTGTCTTGGCCGTGG
AGGTGGCCACCCATTACCTCGAGATCATAACCCAGCTTATAGTGGAGAGCTTCCACTTCAAGAAATGGAGA
AGATGCCCCAGATCTTCTGAAAGTCATCACTAAGCCCTTCAAAAAGCTCATTGTCCAGCTGGATAAAAA
GTTATCAGCCAAATTGCAATGAACGATGAAAAAGCGAAAAACAAGAGTCTTGTCAAGATTGGTGCAAAA
CTTTTACCAACAAGACCCAGATTAACGTCACCTGTTCCCTCGACTGCTAACTGCACCTCCCCTTCCCTCTG
TTGGACGGATGGCATCCAAAACCTGGACCATGAAGAATGTACCTACAAGGAGAACATCGCCAAATGCCAG
CATATCTTTGTGAATTTCCACCTCCCGGATCTTGCTGTGGCACCATCTGTCTATACTCTCCCTGCTGG
TCCTCTGTGGTTGCCTGATCATGATTGTCAAGATCCTGGGCTCTGTGCTCAAGGGGCAGTCCACTGT
CATCAAGAAGACCATCAACACTGATTTCCCTTTCCCTTTGCATGGTTGACTGGCTACCTGGCCATCCTC
GTCGGGGCAGGCATGACCTTCATCGTACAGAGCAGCTCTGTGTTACGTCGGCCTTGACCCCTGATTG
GAATCGGCGTGATAACCATTTGAGAGGGCTTATCCACTCACGCTGGGCTCCAACATCGGCACCACCACCAC
CGCCATCCTGGCCGCTTAGCCAGCCCTGGCAATGCATTGAGGAGTTCACTCCAGATCGCCCTGTGCCAC
TTTTTCTTCAACATCTCCGGCATCTGCTGTGGTACCCGATCCCCTTCACTCGCTGCCATCCGCATGG
CCAAGGGGCTGGCAACATCTCTGCAAGTATCGCTGGTTCGCCGCTTCTACCTGATCATCTTCTTCTT
CCTGATCCGCTGACGGTGTGGCCTCTCGCTGGCCGGCTGGCGGTGCTGGTTGGTGTGGGGTCCC
GTCGTCTTCATCATCATCTGGTACTGTGCTCCGACTCCTGCAGTCTCGCTGCCACGCGTCTGCCGA
AGAACTCCAGAAGTGGAACTTCTGCCGCTGTGGATGCGCTCGTGAAGCCCTGGGATGCGCTCGTCTC
CAAGTTACCCGGCTGCTTCCAGATGCGCTGCTGCTGCTGCTGCCGCTGTGCTGCCGCGGTGCTGCTTG
CTGTGTGGCTGCCCAAGTGTGCGCTGCAGCAAGTGTGCGAGGACTTGGAGGAGCGCAGGAGGGGC
AGGATGTCCCTGTCAAGGCTCCTGAGACCTTTGATAACATAACCATTAGCAGAGAGGCTCAGGGTGAGGT
CCCTGCCTCGGACTCAAAGACCGAATGCACGGCCTTG

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC230431 representing NM_001177999
Red=Cloning site Green=Tags(s)


```
MAPWPPELGAQPNPKYLEGAAGQQPTAPDKSKETNKNNTEAPVTKIELLPSYSTATLIDEPTEVDDPW
LPTLQDSGIKWSERDTKGKILCFQIGRLILLGFLYFFVCSLDILSSAFQLVGGKMQGQFFSNSSIMS
NPLLGLVIGLVTVLVQSSSTSTIVSMVSSLLTVRAAPIIMGANIGTSITNTIVALMQVDRSEFR
RAFAGATVHDFFNWLVLLPVEVATHYLEIITQLIVESFHFKNGEDAPDLLKVIKPFKLIQVLDKK
VISQIAMNDEKAKNLSLVKIWCKTFTNKQINVTVPSTANCTSPSLCWTGDIQNWTKMKNVYKENIAKQ
HIFVNFHLPDLAVGTILLILSLLVLCGCLIMIVKILGSLVKGQVATVIKKTINTDFPFFAWLTGYLAIL
VGAGMTFIVQSSSVFTSALPLIGIGVITIERAYPLTLGSGNIGTTTTAILAALASPGNALRSSLQIALCH
FFFNISGILLWYPIPFTRLPPIRMAGLGNISAKYRWFVAVFYLIIFFFLIPLTVFGLSLAGWRVLVGVGP
VVFIIILVLCRLQLSRCPRVLPKQLQNNWFLPLWMSLKPWDVAVSKFTGCFQMRCCCCCRVCCRACCL
LCGCPKCCRCSCCEDLEEAQEQDVPVKAPETFDNITISREAQGEVPASDSKTECTAL
```

SGPTRRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8073_a07.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001177999

ORF Size: 2067 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001177999.1](#), [NP_001171470.1](#)

RefSeq Size: 4187 bp

RefSeq ORF: 2070 bp

Locus ID: 10568

UniProt ID: [O95436](#)


Cytogenetics: 4p15.2

Protein Families: Transmembrane

MW: 75.6 kDa

Gene Summary: The protein encoded by this gene is a pH-sensitive sodium-dependent phosphate transporter. Phosphate uptake is increased at lower pH. Defects in this gene are a cause of pulmonary alveolar microlithiasis. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, May 2010]

Product images:



Circular map for RC230431