

## Product datasheet for **RC230459L4V**

### PREPL (NM\_001171603) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PREPL (NM_001171603) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PREPL
Synonyms:	CMS22
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001171603
ORF Size:	2181 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC230459).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001171603.1</a> , <a href="#">NP_001165074.1</a>
RefSeq Size:	6318 bp
RefSeq ORF:	2184 bp
Locus ID:	9581
UniProt ID:	<a href="#">Q4J6C6</a>
Cytogenetics:	2p21
Protein Families:	Druggable Genome, Protease
MW:	83.9 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene belongs to the prolyl oligopeptidase subfamily of serine peptidases. Mutations in this gene have been associated with hypotonia-cystinuria syndrome, also known as the 2p21 deletion syndrome. Several alternatively spliced transcript variants encoding either the same or different isoforms have been described for this gene.[provided by RefSeq, Jan 2010]